## MARQUESINA - M.1

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# Desing Features

## 1. - GENERAL RULES:

Being appliances aimed at the elderly group, all elements have been designed following the provisions of Directive 93/42/EEC and RD414/96 Spanish transposition.

Studied and designed by rehabilitation specialists.

Designed to prevent damage and vandalism problems.

Modules accompanied by a panel explaining the correct handling of the appliance, signed by a professional in the field of rehabilitation: Dr. D. Fermin Alvaro Pascual (rehabilitation physician - Collegiate No. 13933, member of the Board FACME - Spanish Federation of Associations Scientific).

## 2. - CERTIFICATES:

- "CE" R.D.414/96

- Authorization No. Madrid DASP233 / M

- License of M<sup>o</sup> Health and Consumer No. 1187-PS

- Quality Certificate ISO 9001 with the following scope: "Design, manufacture and marketing of devices, assistive technologies and features for applications in rehabilitation, special education and geriatrics."

## 3. - FINISH AND SAFETY:

All manufactured products intended to seniors, are manufactured and meet the "CE" RD414/96.

a. Nothing has terminations or sharp or pointed areas.

b. No risk of trapping or crushing.

c. The panels are made of HPL instead of wood to prevent the risk of splintering.

d. All areas are non-slip.

## 4. - MANUFACTURER:

ORTOTECESA, S.L.

C /. Peña de Francia,8

28,00-Arganda del Rey (Madrid)

Tel 918719132

Fax. 918719323

E-mail: [info@ortotecsa.com](mailto:info@ortotecsa.com)

Web: [www.ortotecsa.com](http://www.ortotecsa.com)

## General construction Features

### 1. - METAL STRUCTURES AND PARTS:

#### Series 1 and 2:

(Marquees and modules) Structures Galvanized steel and subsequently coated polyester.  
Stainless steel plate in the ceiling and on the front panel.

Elements of exercises of stainless steel and then coated of polyester.

#### Series 3:

Made of steel, pretreated with an anticorrosive coating "ZN" formulated with high strength epoxy resin and then coated on polyester.

### 2. - ANCHORS:

Galvanized steel.

Vandal-resistant stainless steel.

### 3. - BOARDS (*panels and information points*):

High density laminated panels (HPL).

Resistant to moisture, frost, hail, ultraviolet rays, to acid rain, termites, etc..

High impact strength.

Besides these features binds guarantee that a product is not flammable under approval No - CN267B11CD2000004 class 2 and class CN267B11CD100005 No. 1.

Technical data:

- Thickness tolerance (EN 438-2 ISO 45867/II)
- Density (DIN 23479)
- Thermal conductivity (DIN 52612)
- Dimensional humidity at 20 ° variable (EN 438-2 ISO 4586/II)

Based on tests conducted by TÜV SÜD, this product contains no Polyclinics Aromatic Hydrocarbons (PAH), so it is not toxic. The maintenance is easy to perform, are cleaned with detergent and water.

### 4. - MAINTENANCE:

Maintenance and / or repair to be performed during the warranty period:

Frequency: every six months, preferably in spring and autumn.

Check the operation of the equipment.

Check mounting parts and / or anchor.

Repair and / or replacement of the appropriate element.

Lubricate the bearings and friction areas.

Clean panels (HPL) with soapy water or solvent for paints plastic.

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## **BASIC RECOMMENDATIONS FOR PROPER INSTALATION OF AREAS FOR SENIORS**


1. First request the client all the information and updated documentation on the situation in the area in which you are installing, in relation to existing services (sewerage systems, electrical ...), foundations and generally, any other hidden element that may get in the way at the time of executing the work.
2. If the installation is to be made on parking decks, besides drawings of the existing services, data must be provided on the resistance, depth of the deck, the possibility or not to drill on the same, location of pillars and cross arms ... also because in many cases the use of special anchors is required.
3. In earthen areas, it is recommended to perform a leveling of the land, conferring a slope of approximately 1% to the outside, so as to avoid the accumulation of water in the central areas.
4. Due to the wear that the areas covered with sand suffer for the use, appear holes, gutters and others, which may also favor the accumulation of water. For this reason and for hygiene, they must be raked and refilled at least twice per month, in addition to change the sand at least 1 time per year.
5. It can be installed on the entire surface a solid wood floor framing, with its corresponding treatment, in order that it bears perfectly the weather, with a minimum maintenance. The mentioned floor must be installed so the anchors of the elements remain embedded in the ground, never surfacing over this network to avoid obstacles in the area of use of the exercise elements.
6. Consider leave enough space for the use of each of the elements, according to the destination of the same, also allowing access for people with physical disabilities.
7. The foundation must be realized following the installation plan of each element and so that the concrete never surface to the exterior. For this reason it is recommended that in the earthen area the concrete remains at least 3 centimeters deeper than the land line, as long as the field allows it.
8. The superior part of the concrete cubes must be finished in round to avoid the presence of sharp edges if for any reason outcrops occur. The inferior part of the concrete cube is finished in trapezoidal shape.

### **VERY IMPORTANT:**

**ALL THE ANCHORS THAT FIXED ELEMENTS TO THE GROUND, MUST BE SECURED WITH A FEW SPOT WELDS TO PREVENT VANDALISM.**

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
**STEP 1st.** Pierce on the place to install 1 hole according to the graphic 1600x600x600mm and fill it with concrete without leaving any sharp edges.


**STEP 2nd.** Place the template attached on the concrete base being flushed with the land line.

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
**STEP 3st.** Once the concrete has set, place the marquee on the rods and screw the nuts M16 giving a spot weld in each one once they are tightened.

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OptionB: If the floor where is going to be installed is made of concrete, replace the hooks for 12 metal anchor screws M16x145mm. Spot welding joining nut and concrete screw.

**SE RECOMIENDA NO INSTALAR EN SUELOS ARENOSOS  
IT IS RECOMMENDED NOT TO INSTALL IN SANDY SOIL**