## AUTOPASIVOS

**PARA PUERTA - F.25/3**

**PARA PARED - F.25/4**

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Focused on re-education of the muscles of the upper extremities.


*RECOMMENDATIONS FOR USE*

**Starting position:**

Patient seated on a stool, by matching as much as possible each shoulder with the vertical rope that lowering the device.

Grab the handles with both arms. While one shoulder is at maximum deflection the other is at maximum extension (see figure).

*Figure 110.* Pulleys for deflection-extension of the shoulder. A: at the start of the movement and B: at the end of the movement.

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
**Execution:**

Raise one shoulder as high as possible (no pain). One hand pull the handle while the other arm goes up and vice versa.

Starting from the same position we can do:

- using weights
- or one arm pulls and the other goes up. Relax and return to the start position.

*Figure 70.* Position for deflection-extension exercises of the shoulder.

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Figure 71. Position for shoulder exercises of deflection-abduction, extension-abduction

Figure 72. position for deflection-internal rotation and extension-external rotation exercises in sitting shoulder.

NOTE: Follow carefully all instructions given by your doctor, physic or therapist for the use of the product.

#### RECOMMENDATIONS FOR INSTALLATION