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ANKLE MOBILITY BASE – F.40/42

Mechanised adjustable base, especially recommended for re-educating the ankle.

Allowing for plantar flexion, dorsiflexion, internal and external rotation.

For resistance exercises, the following weights are provided:

Three x 2 kg., one x 1 kg. and one x 1/2 kg.

Metal frame with chrome and epoxy finish.

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OPTIONAL ACCESORIE:

SEAT FOR ANKLE BASE - F.40/42-1

Adjustable height.
Equipped with weights stand.
Metal frame with chrome and epoxy finish.
Fireproof upholstery.

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RECOMMENDATIONS OF USE

Apparatus for ankle exercise rehabilitation in their variations of plantar flexion, dorsiflexion, internal and external rotation.

The exercises must be ruled by the physiotherapist and be adapted by pathology and evolutionary progress of the patient.

The exercises must be done slowly and in a controlled way.

Exercises for dorsiflexion and plantar flexion

Sitting in a chair that allows us to have the knee in 90 degrees of flexion, put your foot on the center and will make movements into flexion and extension (forward and backward) throughout the possible range, trying not to move the knee as far as possible. we will also make the movement sideways through the movement of the ankle.

Exercise for lateral pronation and supination

Starting from the basic position of the previous exercise, try to bring the sides of the bowl almost to the ground (as in the example of the photo). For that it is important to do with a movement of the ankle and not through the knee, which we will try to maintain fixed.

Exercise for circumduction

We will make circles with the feet from one side to another.

NOTE: if any of the exercises increase or produce excessive pain or swelling, should be stopped and ask the therapist.

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