
# BAÑO DE PARAFINA (22L.) G.18907

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Note:  $U_t$  is the AC mains voltage prior to application of the test level.

<b>GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY GUÍA Y DECLARACIÓN DE ORTOTECSA – INMUNIDAD ELECTROMAGNÉTICA</b>			
The paraffin bath is intended for use in an electromagnetic environment specified below. The customer or user of the paraffin bath should assure that it is used in such an environment			
<b>Immunity Test Ensayo de inmunidad</b>	<b>IEC 60601 test Level Nivel de ensayo de la norma</b>	<b>Compliance Level Nivel de conformidad</b>	<b>Electromagnetic Environment Guidance Entorno electromagnético, guía</b>
<p>Conducted RF IEC 61000-4-6 RF conducida</p> <p>Radiated RF IEC 61000-4-3 RF radiada</p>	<p><math>3 V_{rms}</math> 150 kHz to 80 MHz</p> <p><math>3 V/m</math> 80 MHz to 2,5 GHz</p>	<p><math>3 V_{rms}</math></p> <p><math>3 V/m</math> 80 MHz to 1 GHz (due to the testing laboratory limitations) (debido a las limitaciones de laboratorio)</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the infrared lamp, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Los equipos móviles y comunicaciones de RF no se deben usar más cerca de cualquier parte de la lámpara, incluyendo los cables, que la distancia de separación recomendada a la frecuencia del transmisor.</p> <p>Recommended separation distance: Distancia de separación recomendada:</p> <p><math>d = 1,17 \sqrt{P}</math> 150 kHz to 80 MHz <math>d = 1,17 \sqrt{P}</math> 80 MHz to 800 MHz <math>d = 2,33 \sqrt{P}</math> 800 MHz to 2,5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Donde <math>p</math> es la maxima potencia de salida asignada de transmisor en vatios (w) conforme al fabricante del transmisor y <math>d</math> es la distancia de separación recomendada en metros.</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range. <sup>b</sup></p> <p>Las intensidades de campo desde el transmisor fijo de RF, según se determina por un estudio electromagnético de lugar, debe ser menor que el nivel de conformidad en cada rango de frecuencia.</p> <p>Interference may occur in the vicinity of the equipment marked with the following symbol: La interferencia puede ocurrir en la vecindad del equipo marcado con el siguiente símbolo:</p> 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the			

electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the infrared lamp is used exceeds the applicable RF compliance level above, the XXXX should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation the infrared lamp.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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The paraffin bath is intended for use in a electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the bath can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitter) and the paraffin bath as recommended below, according to the maximum output power of the communication equipment.

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For transmitters rated at a maximum output power nor listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 to 80 MHz and 800 MHz, the separation distance is applied in the highest frequency range.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people.

# PARAFINE BATH (22L.)

## G.18907

### Revision History:

Rev	Reason	Date
1	Creation of the document	17- jun - 2013
2	The inclusion of comments by electrical safety tests	05 -ago - 2013
3	The inclusion of the reference product that in 115Vac	28 -Abr- 2014
4	TF update review by SGS	06- Sep - 2016
5	Drainage recommendation	04- May - 2018

<b>INDEX</b>	<b>PAGE</b>
MODELS	19
EQUIPMENT CHARACTERISTICS	19
ELECTRICAL SPECIFICATIONS	19
WARNING	21
INDICATIONS OF PARAFFIN	22
CONTRAINdicATIONS AND PRECAUTIONS	22
ELEMENTS	24
INSTRUCTION FOR USE	24
APLICACION	25
PRODUCT STORING	25
CLEANING PRODUCT	26
TROUBLESHOOTING GUIDE	26
LABEL ON PRODUCT	27
WARRANTY	28
EMC COMPATIBILITY TABLES	29
MODELS	19

## Models

<u>REF</u>	<u>NAME</u>
G.18907	Paraffine bath for Mains Input Voltage of 230Vac

## Equipment characteristics

Paraffin bath with heating system based on "Bain Marie". For hospital use.

Mobile (lock wheels).

Built both the inner tank and the outer bowl in stainless steel AISI 316.

Edge insulation (wood).

Safety mechanism against boiling based on Thermostat.

Thermostat from 0 to 60 °. Thermostatic control paraffin maintains the optimum temperature for maximum therapeutic benefit.

Front panel provided with pilot and thermostat control.

Optional accessories:

- G.18907 - 1: paraffin bath trolley
- G.18907 - 2: filter for paraffin
- G.18910: 1kg paraffin tablet

## Electrical specifications

### Rating

Models		G.18907
Voltage	110 Vac ≈	230 Vac ≈
Frequency	60 Hz	50 Hz
Power	1800 W	1800 W
Amperage	18A	9,1A
Useful life	10 Years	10 Years

### Main cordon:

- 4 m longitude (13 ft)
- Cordon and plug with ground connection.

### Leaking current

- Conform to IEC 60601-1

## Environment Conditions

### Functional environments conditions recommended

- Temperature: 10 °C - 35 °C (Temperature recommended: 20 °C - 25 °C)
- Humidity: 35 - 75%
- Pressure: 1000 - 700 mBar

**Storage and transport Conditions**

- Temperature: -5 °C - 70 °C
- Humidity: 0% - 90%
- Pressure: 1000 - 100 mBar

**Environment requirements**

- according to directive 2011/65/UE (RoHS2).
- according to directive 2002/96/EC. (WEE).

**WARNING!**

Read the instructions carefully before handling the product.

- Before connecting the product, check that your mains voltage corresponds to that is indicated on the label.
- Do not use the equipment if it does not work properly or is damaged, or if the power cable or plug is defective.
- Never immerse the product in water or use in close proximity to it.
- Do not use the equipment near a bathtub, shower or with wet hands. Should any part of the product was wet, reject the use.
- For cleaning proceed as it is mentioned in conservation section of this manual.
- Do not cover the equipment during use.
- Dont stack bath with other equipment.
- Use the bath on flat and stable, preventing the unit from falling.
- Keep a safe distance bigger than 80 centimeters, between baths and potentially flammable objects.
- Turn off the bath and unplug from the mains before cleaning or maintenance.
- No treatment must applied after taking analgesics or painkillers, ther heat sensitivity may be diminished by the effect of drugs and this can cause skin lesions.
- The ingestión of alcohol or drugs can also alter the sensitivity to heat.
- The device is not intended for persons insensitive to heat.
- Excessive exposure to heat can cause burns or skin lesions.
- Do not use the equipment in areas of the body with burns, infections, inflammations, rashes, wounds or hypersensitive areas.
- Do not sleep with the product connected.
- Persons with disabilities, disabled, handicapped, insensitive to heat, or sick children should not use this product without direct adult supervision.
- It should not be used under the influence of alcohol, sleeping pills or sedatives.
- If you are carrying a pacemaker or similar cardiovascular device, we recommend you always check with your doctor before using the product.
- Should you experience pain or unpleasant sensation during use, stop using and consult your doctor.

**INDICATIONS OF PARAFFIN:**

The wax used in therapy must be white, odorless, tasteless and solid and is supplied in the form of plates.

Therapeutic effects are found with the application of heat:

- To Increase connective tissue extensibility.
- To decrease joint stiffness.
- Analgesic effects
- In subacute and chronic inflammatory processes. Example, rheumatoid arthritis, bursitis, tendinitis and peri-arthritis escapulahumeral, because it provides a long anti-inflammatory and analgesic action.
- Muscle relaxant effect.
- In painful processes as sacrolumbalgias, back pain, cervical pain, sciatic pain, herniated discs, etc..
- Antispasmodic effects.
- Anti-inflammatory effect.
- The efficacy of thermotherapy applications is clinically documented to fight against muscle spasms of defense, subacute and chronic traumatic events and other processes like sympathetic dystrophy - reflects, Raynaud's disease, pelvic inflammatory disease, etc..

**CONTRAINDICATIONS AND PRECAUTIONS****Disorders in the sensitivity.**

Before applying, you must evaluate the thermal and pain sensitivity of the area as well as the patient's circulatory status, because the heat is generally contraindicated or should be used with particular caution on areas anesthetized or obtunded patient.

In most applications, dosimetry is not exact, and must be taken of the patient perceives discomfort or pain when warming threshold was exceeded.

**Severe circulatory disorders.**

Heating of tissues with inadequate blood supply is contraindicated, as the temperature elevation increases the metabolic demands, without associated vascular adaptations, with the risk of producing an ischemic necrosis.

Circulatory disorders, especially arterial, can provide to a decrease in convective mechanism supplied heat dissipation, with the risk of burning the tissues of the area. Altered vessels may have an increased contractility and their response to intense thermal stimulus is essentially spasmodic.

**Bleeding disorders.**

Don't treat areas where there is a tendency to bleed, because the thermal stimulus produces an increase in blood flow. This is the case in diseases such as hemophilia, in post-traumatic states or in patients undergoing long-term steroid therapy, in which there is a tendency to capillary fragility.

**Malignant neoplas.**

Except for special treatment (hyperthermia anticancer), thermotherapy is contraindicated in areas where malignancy is localized, as subtherapeutic temperatures, as used in physiotherapy, can increase the rate of tumor growth. There are many studies showing that thermotherapy applications may increase the development of metastases.

**Pregnant.**

**Febrile states.**


Avoid applications over recent injuries, whether or not infected by the production risk of burns. Be very careful in the application of recent and skin grafts, and, similarly, on infected skin lesions (risk of exacerbation of these) in peripheral vascular disease or any condition that affects the distribution of internal heat.

It is important to supervise  
the level of the water and that  
it is the proper one.

**ELEMENTS:**

1. Tank of Paraffin
2. Thermostate
3. Pilot
4. Switch to ON/OFF
5. Water level indicator
6. Wheel brake
7. Mains cable connetion
8. Hole to fill the water tank
9. Cap water tank
10. Lid of paraffin tank
11. paraffin tank drain
12. drain water tank
13. Reset button

**INSTRUCTION FOR USE:**

Before connecting the unit to the network is necessary to fill the tank with water.

Fill the water tank (8). Check that the tank is full to capacity by the water level indicator (5), it must be filled completely.

With the water tank full of water, fill the tank (1) with paraffin, insert the plug (7) to the Mains and turn the switch (4) to the position 1 (On).

The pilot light (3) will ON.

Select the temperature for the paraffin desired for the thermostat (2).

When the pilot turns OFF, it had reached the temperature selected.


We recommend a 50-55°C temperature.

The warning pilot, will blink from time to time, indicating that the resistance that heats the water in the tank is working to maintain the temperature selected on the thermostat.

Both models of bath are equipped with a system that cuts power before the water reaches 100 ° C. in case of failure of the main thermostat.

If Mains accidentally was cut, the unit can be reset after a period of about 15 minutes for it to access the bottom of the unit, and press the reset button (13).

If a sign of failure arises again, in the main thermostat, so you should call Customer Service.

**Attention!**

Do not stack to another equipment or device.

Using distilled water for good maintenance of the paraffin bath is recommended.

To proceed to the emptying tank, access outfalls paraffin (11) and water (12).


**Caution:** Periodically check the water level. The implementation team with insufficient water level can damage the heating elements.

SELECTOR FOR TEMPERATURE  
OF PARAFFINE. Thermostat.

## CONTROL PANEL

**Important Notice!:** The temperature of the wax will be selected by the thermostat (2), does not correspond with the temperature dial indicator water level (5), since it does not act as a thermometer (this will mark a difference to least about 20 °).

## APPLICATION

- Before treatment, the body part should be cleaned with soap and water, then with alcohol to remove any soap residue and prevent bacterial growth on the bottom of the bath.
- Carefully pour the wax with a casserole on the hand or foot or put your hand or foot for several seconds in the bathroom. Then take off, in order to form a thin coating of paraffin slightly hardened and adherent to the skin .
- The operation is repeated 8 to 12 times until it forms a thick layer of paraffin wax. Then the area is wrapped in a plastic bag and covered with several towels in order to facilitate heat retention.
- The body area should be removed of any metallic objects and must ensure that the area does not move once introduced in the bath to avoid the appearance of "hot spots".
- By moving the part introduced in the molten paraffin, there is a risk of interrupting the semisolid paraffine barrier, by which the patient could feel a burning sensation. The hand will be immersed with fingers extended and separated as possible.
- The patient should be in a comfortable position, with elevated until the end of the treatment, to avoid the potential appearance of edema.
- The application is maintained for 15 to 20 min.
- After this time, remove the towel and bag, and remove the layer of solid paraffin and put into the bath, but previously sterilized at 80 ° C, or with boiling water
- After application, verify the condition of the skin.
- After application, the area should be cleaned with soap and water. The cleaning can be completed with a gentle massage with a moisturizing lotion or mineral oil to moisturize and soften the skin.
- After applying wax, the skin is smoother, softer, moist and flexible, making it easier to massage and mobilization.

## PRODUCT STORING

- To extend the life of the product, you must unplug it and keep it disconnected from the outlet when not be used for some periods of time.
- Store in a dry place.
- Allow to cool before storing.
- Don stack objects in its top.
- Do not place heavy objects.

**CLEANING PRODUCT**

The bath needs continued supervision to avoid contamination.

The bath should be cleaned regularly, as recommended by ORTOTECSA. The times you reuse the paraffin determine the frequency with which the bath is to be cleaned, but it is recommended in intervals not exceeding six months.

**MAINTENANCE GUIDE**

1. Inspect the environmental conditions in which the bath is located. Avoid storage in wet places. Do not stack other devices or equipment on top of the paraffin bath.
2. Perform integral external and internal cleaning service.
3. Externally inspect the equipment and its main accessories.
4. Check cord and other electrical components.
5. Lubricate moving parts with any commercial lubricant spray.
6. Check the accuracy of temperature control.
7. Check the correct operation of the equipment in conjunction with the operator.
8. The glove or sock which is formed during the treatment should not be added directly to the bath, but collect and sterilized at 80 ° C, or with boiling water.
9. It is also necessary to renew entirely contents of the bath every six months and clean it properly, since the surface with traces of skin is an excellent medium for bacterial culture.
10. In physiotherapy services it is recommended to have at least one unit for hands and one for feet.
11. The unavoidable dirt residues originate by the paraffin in the around zones and in table, on the floor etc., should put the units in a separate area.

**PARAFFIN CLEANING SYSTEM**

- Paraffin cleaning is performed according to their use and degree of dirtiness.
- The temperature is increased to proceed to the emptying of the bath (Recommended 90°).
- Empty the bath content in a tank for this type of residues according to the laws in force in the country. Don't pour down into the drain or down the toilet.
- The procedure to clean the bathroom with soap and water. It should be necessary to pay attention in all those areas with corners and more difficult to reach. If it is necessary use a brush.

**TROUBLESHOOTING GUIDE:**

Possible malfunctions.

INCIDENCE	RECOMMENDATION
Bath doesn't switching-on	<p>Check that there is right voltage at the socket used. Press the switch to switch-ON .</p> <p>Check that the thermal cut-out has not tripped (this can jump caused by paraffin overheating or caused for not maintain the water level of the tank at the level indicated in the instructions).</p> <p>Check that the equipment is used with the proper water level (otherwise, it can tripped the thermal cut-out and burn resistance). Check the conditions of the power cord.</p> <p>If still remains the problem go to our authorized service.</p>
Bath doesn't heat.	Si el baño de parafina está encendido pero no

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**MEANING OF SYMBOLS**

SYMBOL	DESCRIPTION
	MANUFACTURER
	MANUFACTURE DATE (YEAR/MONTH)
<b>MODELO:</b>	REFERENCE MODEL OF DEVICE
<b>SN:</b>	SERIAL NUMBER
<b>LOTE:</b>	PART NUMBER OF DEVICE
	CONSULT THE ACCOMPANIMENT DOCUMENTATION
	The device must not be disposed as an urban waste, because it is an Electrical and electronic waste
	TECHNICAL DATA

<b>CE</b>	CE marking
-----------	------------

**WARRANTY**

ORTOTECSA, S.L. will provide a warranty for two years from the date of purchase against manufacturing defects. Our warranty includes repair, replacement, or exchange of the product and / or components including labor in our installations.

To have effect this guarantee, is required the purchase invoice.

To have effect this guarantee, is required the purchase invoice.

**THIS WARRANTY SHALL NOT UNDER THE FOLLOWING CONDITIONS:**

- When the use, care and operation of the product has not been done in accordance with the instructions contained in the operating instructions.
- When the product has been used outside of its capacity, abused, beaten, exposed to moisture, ingressed by any liquid or corrosive substance, as well as any other fault attributable to the consumer.
- When the product has been disassembled, modified or repaired by unauthorized persons.
- The breakage and / or replacement of the bulb are not part of the warranty for this product.

In our continuous improvement process, Ortotecsa reserves the right to modify any characteristic without prior notice.

**EMC EN 60601-1-2-2007****EMISSIONS**

The infrared lamp uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

The infrared lamp is suitable for use in all establishments, including domestic establishments and those directly connected to the mains supply low voltage which supplies buildings used for domestic purposes. It doesn't built any active electronic component.

**INMUNITY**

Floors should be of wood, concrete or ceramic. If floors are covered with synthetic material, the relative humidity should be at least 30%.

The unit is tested to be immune to electric fields of:

- $\pm 6$  kV by contact.
- $\pm 8$  kV by air.

The quality of the Mains should be that of a typical commercial or hospital.

Normally lamps will support without show any problema of working:

- Quick bursts:  $\pm 2$  kV for Mains.
- Shockwaves:  $\pm 1$  kV differential mode

The quality of the power supply must be that of a typical commercial or a hospital. If the user of the infrared lamp requires continued operation during power interruptions, it is recommended that the lamp be powered from an uninterruptible power supply or a battery.

Magnetic fields at Mains frequency should be like a typical location in a commercial environment or hospital: 3 A / m.

Mobile equipment and RF communications should not be used closer to any part of the lamp, including cables. The recommended separation distance to the transmitter frequency is depending of the working freq of the device:

- $d = 1.17 \sqrt{P}$  150 kHz to 80 MHz
- $d = 1.17 \sqrt{P}$  80 MHz to 800 MHz
- $d = 2.33 \sqrt{P}$  800 MHz to 2.5 GHz

Where  $P$  is the maximum output power of the transmitter in watts (w) according to the transmitter manufacturer specification and  $d$  is the recommended separation distance in meters.

Field strengths from fixed RF transmitter, as determined by an electromagnetic site inspection and should be less than the compliance level in each frequency range of the above table.

Interference may occur in the vicinity of equipment marked with the following symbol:

See tables attached to the end of the manual.

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			vecindad del equipo marcado con el siguiente símbolo: 
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies          NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the infrared lamp is used exceeds the applicable RF compliance level above, the <b>XXXX</b> should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation the infrared lamp.  <sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

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The parafine bath is intended for use in a electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the infrared lamp can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitter) and the infrared lamp as recommended below, according to the maximum output power of the communication equipment.

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For transmitters rated at a maximum output power nor listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 to 80 MHz and 800 MHz, the separation distance is applied in the highest frequency range.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people.

# ORTOTECSA,S.L.

Fabricante/Manufacturer: Ortotecsa,s.l.  
C/Peña de Francia,8 – 28500 Arganda del Rey (Madrid), España.  
Tlf.: +34918719132 / Fax: +34918719323  
email: [info@ortotecsa.com](mailto:info@ortotecsa.com), [www.ortotecsa.com](http://www.ortotecsa.com)