CALENTADOR COMPRESAS

MEDIANO – G.18912/M

GRANDE – G.18912/L

EXTRA GRANDE - G.18912/XL

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LARGE – G.18912/L

EXTRA LARGE – G.18912/XL

Made from stainless steel AISI 316 both interior and exterior.

Thermostat ranging to +130°C.

Recommended temperature: 60°C.

3 figure viewing display.

All the work parameters can be adjusted using the buttons at the front and are stored in the memory.

Mobile /equipped with handles and casters).

Includes a stainless steel draining rack with a capacity to hold 6-6 compresses.

Degree of protection: IP20 for the whole apparatus, IP52 for the front part.

If the thermostatic probe brakes, the heater will automatically be switched OFF.

Electric characteristics: 230V/50 Hz.

Electric characteristics: 120V/60 Hz. (Models 120V)


This machine compilies with health product standard MDD 93/42 CEE.

INSTRUCTIONS

Fabricante: Ortotecsa,s.l. C/Peña de Francia,8 – 28500 Arganda del Rey (Madrid), España.
Tlf.: +34918719132 / Fax: +34918719323 – email: info@ortotecsa.com, www.ortotecsa.com

COMPONENTS

- 1 / START SWITCH
- 2 / VISOR WATER TEMPERATURE AND PROGRAMMER
- 3 / SUPPORT FOR COMPRESS
- 4 / BUCKET
- 5 / ELECTRICAL RESISTANCE
- 7 / CABLE CONNECTION TO MAINS
- 8 / KEY FOR DRAIN (ON A SIDE)

INSTALLATION AND COMMISSIONING

IS SUITABLE INSTALL THE HEATER WHERE THERE IS A DRAIN NEAR.
TO RENEW THE WATER OPEN THE COCK (8).


BEFORE RUN THE HEATER, FILLING UP THE BOWL (4) WITH 10 GALLON WATER.
YOU SHOULD WATCH THE WATER LEVEL, TO DON'T LOSE 60mm, COUNTING FROM THE
BOTTOM OF THE BUCKET, SO THAT COVERS THE RESISTANCE (5) AND REMAIN BELOW
THE SUPPORT FOR COMPRESS.

TEMPERATURE RECOMENDED 60°C.

DESCRIPTION OF THE FRONTAL PANEL (32x75mm):

- 1 – Key Set: allows programming the Set Point and the operating parameters
- 2 - Key **DOWN** : allows decreasing the preset values and the selection of the parameters
- 3 - Key **UP/DEFROST** : allows increasing the preset values, to select the parameters and to activate the defrost manual
- 4 - **Key F**: It can be programmed using the parameter “USrb to actuate as ON/OFF (Stand-by). In the programming mode of the hidden parameters, it can be used to modify the visibility of the parameters (see PROGRAMMING MANUAL)
- 5 - **Led SET**: Led SET: indicates that we are in the programming mode and the programming level of the parameters. It also indicates that the equipment is in Stand-by mode.
- 6 - **Led OUT**: indicates the state of the output of the compressor (or the temperature control device): ON, OFF or INHIBITED (flashing)

CONTROL TEMPERATURE

Programming the required temperature (Set)

1. Press “SET” and release. The display will show SP (Set Point)

2. Press “UP” or “DOWN” to reach the required value. If you hold the key UP or DOWN the speed increase quickly

3. Press “SET” to confirm the value

On and OFF by keypad

Press “F” for 2 seconds. To activate this function, program the parameter “uSrb=1” (see access to hidden parameters)

ACCESS TO STANDARD PARAMETERS

Press “SET” for 5 seconds. Display will show the first visible parameter (usually SPLL) (if shows 0, insert password number and press “Set”)

2. Press “UP” or “DOWN” to reach the required parameter

Programming the parameter

1. Press "SET" to enter the parameter.

2. Press "UP" or "DOWN" to change the value of the parameter

3. Press "SET" to confirm the value

INDEX OF STANDARD PARAMETERS

Rango		Def.	
SLS	Lower bound of the Set Point	-58.0÷SPHL	
ISE	probe configuration	PTC-NTC	PT
IUP	Decimal point	C0,F0,C1,F1	C1
r.d	Differential (hysteresis)	0.0÷30.00 °C/°F	2.0
r.Hc	Relay operating model	Heat-Cool	Heat

TEMPERATURE REGULATION

The thermostat setting mode is ON / OFF and acts on the output OUT (relay) based on the reading of the probe, on the Set Point "SP" (value of the required temperature), on the differential "HSEt (hysteresis) and on the operating mode "Func".

Depending on the operating mode set in the parameter "Func", the differential is considered automatically positive to control for cold (Func = Cool) and negative to control for heat (Func = HEAt); eg chilled up to 5 degrees SP and return to start the cold to 8 degrees SP + HSEt.

In case of probe error you can make the output OUT continue to work in cycles as scheduled in the "tonE" parameter (output activation in minutes) and "toFE" (strike out in minutes).

In probe error condition, if "toFE=OFF" and "tone= to a value" then the output OUT will remain activate all the time

INDEX ALARM AND SIGNAL

E	failure, configuration parameter SEnS, damage or out of range. In case of probe error, the controller works for the times programmed in the parameters tonE and toFE
EEPR	Memory error (turn ON and OFF the equipment)
Od	Delay relay (as scheduled in the parameter Od)
dEF	Desfrost activate
PdEF	Final defrost to return to the set temperature (SP+Etdv)

ELECTRICAL SCHEMATIC