BAÑOS DE REMOLINOS

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última revision 07/08/2019

BAÑOS DE REMOLINOS	BAÑOS ESTANDAR	
	G.18951/A	MOVIL 60L (15 GAL) - 220V
	G.18951/A-1	ALTURA AJUSTABLE 60L (15 GAL) - 220V
	G.18951/B	MOVIL 80L (22 GAL) - 220V
	G.18951/B-1	ALTURA AJUSTABLE 80L (22 GAL) - 220V
	G.18951/C	MOVIL 100L (27 GAL) - 220V
	G.18951/D	MOVIL 170L (45 GAL) - 220V
	G.18951/E	MOVIL 320L (85 GAL) - 220V
	G.18951/X	BAÑOS A MEDIDA
	G.18951/40	MOVIL 40L (10 GAL) - 220V
	G.18951/230	MOVIL 230L (60 GAL) - 220V
	G.18951/290	MOVIL 290L (75 GAL) - 220V
	G.18951/340	MOVIL 340L (90 GAL) - 220V
	G.18951/380	MOVIL 380L (100 GAL) - 220V
	G.18951/400	MOVIL 400L (105 GAL) - 220V
	G.18951/1000	ESTACIONARIO 1000L (260 GAL) – 220V
G.18951/1140	ESTACIONARIO 1140L (300 GAL) – 220V	

BAÑOS DE REMOLINOS 120V	BAÑOS ESTANDAR	
	G.18951/A (120V)	MOVIL 60L (15 GAL) - 120V
	G.18951/A-1 (120V)	ALTURA AJUSTABLE 60L (15 GAL) - 120V
	G.18951/B (120V)	MOVIL 80L (22 GAL) - 120V
	G.18951/B-1 (120V)	ALTURA AJUSTABLE 80L (22 GAL) - 120V
	G.18951/C (120V)	MOVIL 100L (27 GAL) - 120V
	G.18951/D (120V)	MOVIL 170L (45 GAL) - 120V
	G.18951/E (120V)	MOVIL 320L (85 GAL) - 120V
	G.18951/X (120V)	BAÑOS A MEDIDA
	G.18951/40 (120V)	MOVIL 40L (10 GAL) - 120V
	G.18951/230 (120V)	MOVIL 230L (60 GAL) - 120V
	G.18951/290 (120V)	MOVIL 290L (75 GAL) - 120V
	G.18951/340 (120V)	MOVIL 340L (90 GAL) - 120V
	G.18951/380 (120V)	MOVIL 380L (100 GAL) - 120V
	G.18951/400 (120V)	MOVIL 400L (105 GAL) - 120V
	G.18951/1000 (120V)	ESTACIONARIO 1000L (260 GAL) – 120V
G.18951/1140 (120V)	ESTACIONARIO 1140L (300 GAL) – 120V	

Contenidos

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			<p>conformidad en cada rango de frecuencia. Interference may occur in the vicinity of the equipment marked with the following symbol: La interferencia puede ocurrir en la vecindad del equipo marcado con el siguiente símbolo:</p> 
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation the unit.</p>			
<p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

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NOTE 1 to 80 MHz and 800 MHz, the separation distance is applied in the highest frequency range. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people.

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Last revision 22/09/2021

Contents

1.	Introduction	28
1.1.	About this manual	28
1.2.	Important information	28
2.	Manufacturer	29
3.	Purpose of product	29
4.	Standar Models and general dimensions.....	30
4.1.	Models working in 220VAC	30
4.2.	Draw and dimensions models working in 220VAC.....	30
4.3.	Models working in 120VAC	33
4.4.	Draw and dimensions models working in 120VAC.....	33
5.	Classification according to Medical Devices Directive 93/42/EEC and normative EN 60601-1:2006	33
6.	Operator skills and training	34
7.	Patient evaluation	35
8.	Inspecting the whirlpool bath	35
9.	While using the whirlpool	36
10.	Instructions for use	36
11.	Cleaning and disinfecting instructions	38
12.	Maintenance and troubleshooting	39
12.1.	Maintenance Schedule	39
12.2.	Maintenance procedures	39
12.3.	Troubleshooting	39
12.4.	Labelling	40
13.	Information related to EMC	41
14.	Warranty	46

WHIRLPOOL BATHS	BAÑOS ESTANDAR	
	G.18951/A	MOVIL 60L (15 GAL) - 220V
	G.18951/A-1	ADJUSTABLE HEIGHT 60L (15 GAL) - 220V
	G.18951/B	MOVIL 80L (22 GAL) - 220V
	G.18951/B-1	ADJUSTABLE HEIGHT 80L (22 GAL) - 220V
	G.18951/C	MOVIL 100L (27 GAL) - 220V
	G.18951/D	MOVIL 170L (45 GAL) - 220V
	G.18951/E	MOVIL 320L (85 GAL) - 220V
	G.18951/X	BAÑOS A MEDIDA
	G.18951/40	MOVIL 40L (10 GAL) - 220V
	G.18951/230	MOVIL 230L (60 GAL) - 220V
	G.18951/290	MOVIL 290L (75 GAL) - 220V
	G.18951/340	MOVIL 340L (90 GAL) - 220V
	G.18951/380	MOVIL 380L (100 GAL) - 220V
	G.18951/400	MOVIL 400L (105 GAL) - 220V
	G.18951/1000	STATIONARY 1000L (260 GAL) – 220V
G.18951/1140	STATIONARY 1140L (300 GAL) – 220V	

WHIRLPOOL BATHS 120V	BAÑOS ESTANDAR	
	G.18951/A (120V)	MOVIL 60L (15 GAL) - 120V
	G.18951/A-1 (120V)	ADJUSTABLE HEIGHT 60L (15 GAL) - 120V
	G.18951/B (120V)	MOVIL 80L (22 GAL) - 120V
	G.18951/B-1 (120V)	ADJUSTABLE HEIGHT 80L (22 GAL) - 120V
	G.18951/C (120V)	MOVIL 100L (27 GAL) - 120V
	G.18951/D (120V)	MOVIL 170L (45 GAL) - 120V
	G.18951/E (120V)	MOVIL 320L (85 GAL) - 120V
	G.18951/X (120V)	BAÑOS A MEDIDA
	G.18951/40 (120V)	MOVIL 40L (10 GAL) - 120V
	G.18951/230 (120V)	MOVIL 230L (60 GAL) - 120V
	G.18951/290 (120V)	MOVIL 290L (75 GAL) - 120V
	G.18951/340 (120V)	MOVIL 340L (90 GAL) - 120V
	G.18951/380 (120V)	MOVIL 380L (100 GAL) - 120V
	G.18951/400 (120V)	MOVIL 400L (105 GAL) - 120V
	G.18951/1000 (120V)	STATIONARY 1000L (260 GAL) – 120V
G.18951/1140 (120V)	STATIONARY 1140L (300 GAL) – 120V	

1. Introduction

1.1. About this manual

This manual is a combined user/operation and maintenance manual. It is applicable to all models found in 4.1 and 4.3

1.2. Important information

Read Working Instructions

Read these instructions carefully before handling the product.

- ▶ If these instructions are not kept, could result in damage to your new whirlpool and/or body injury.
- ▶ Before connecting the product, check that your mains voltage corresponds to that is indicated on the label.
- ▶ Do not use the equipment if it does not work properly or is damaged, or if the power cable or plug is defective
- ▶ Whirlpool baths are intended just for professional use only.
- ▶ DO NOT operate appliance without properly filling with water.
- ▶ Fill appliance with water to the prescribed level before plugging into an electrical receptacle.
- ▶ Operating the appliance without water may result in damage to the motor.
- ▶ Improper use of the whirlpool can cause injury. Use the whirlpool only for the purpose described in this manual.
- ▶ The turbine is top-heavy. Improper handling can cause injury or damage. Handle the turbine with care.
- ▶ Connection improper items to the whirlpool baths can cause injury to persons and damage to equipment. Use only Ortotecs approved items in whirlpools.
- ▶ Never immerse the product in water or use in close proximity to it.
- ▶ Do not use the equipment near a bathtub, shower or with wet hands. Should any part of the product was wet, reject the use.
- ▶ For cleaning proceed as it is mentioned in conservation section of this manual.
- ▶ Do not cover the equipment during use.
- ▶ Don't stack bath with other equipment.
- ▶ Use the bath on flat and stable, preventing the unit from falling.
- ▶ Keep a safe distance bigger than 80 centimetres, between baths and potentially flammable objects.
- ▶ Turn off the bath and unplug from the mains before cleaning or maintenance.
- ▶ Do not sleep with the product connected.
- ▶ Persons with disabilities, disabled, handicapped or sick children should not use this product without direct adult supervision.
- ▶ It should not be used under the influence of alcohol, sleeping pills or sedatives.
- ▶ If you are carrying a pacemaker or similar cardiovascular device, we recommend you always check with your doctor before using the product.
- ▶ Should you experience pain or unpleasant sensation during use, stop using and consult your doctor.
- ▶ Don't locate the unit for the treatment in such way be difficult unplug it from the outlet if was necessary by fails
- ▶ When the unit has reached its life time should be disposed according to the local government laws. Please get in contact with the local authorities for more information.

- ▶ There are no special environmental conditions for the use. An indoor environment is perfect for the use. Nevertheless the recommendation is: temperature of 25-28º, humidity between 30-50%, Atmospheric pressure: 1000-1025hPa-mb
- ▶ Don't manipulate or open the unit: it is not allowed except by the manufacturer personnel.
- ▶ The Mains cable must NOT be replaced. Only manufacturer is allowed to do it.

2. Manufacturer

Name: ORTOTECSA S.L
Address: C/Peña de Francia,8
28500 Arganda del Rey
Madrid. España.
Tlf.: +34 918 719 132
Fax: +34 918 719 323
email: info@ortotecsca.com,
WEB: www.ortotecsca.com.

3. Purpose of product

The Whirlpool, called whirlpool in the Anglo-Saxon terminology, consisting of baths whose water is kept in constant agitation by a turbine. In these baths, the thermal effects of hot or cold water are added to derivatives from agitation of stirring water. It is one of the hydrotherapeutic methods most studied and currently used in the treatment of physical dysfunction.

The effects of stirring water bath on tissues have been studied by several authors. The turmoil created in the whirlpool bath works as a source of stimulation mechanical skin, acting as a anti-irritant and encouragement of large sensitive afferents, blocking the transmission of pain. Moreover, the increased agitation convective will propagate of hot mechanism.


4. Models and general dimensions

4.1. Models working in 220VAC

WHIRLPOOL BATHS	BAÑOS ESTANDAR	
	G.18951/A	MOVIL 60L (15 GAL) - 220V
	G.18951/A-1	ADJUSTABLE HEIGHT 60L (15 GAL) - 220V
	G.18951/B	MOVIL 80L (22 GAL) - 220V
	G.18951/B-1	ADJUSTABLE HEIGHT 80L (22 GAL) - 220V
	G.18951/C	MOVIL 100L (27 GAL) - 220V
	G.18951/D	MOVIL 170L (45 GAL) - 220V
	G.18951/E	MOVIL 320L (85 GAL) - 220V
	BAÑOS A MEDIDA	
G.18951/X		
G.18951/40	MOVIL 40L (10 GAL) - 220V	
G.18951/230	MOVIL 230L (60 GAL) - 220V	
G.18951/290	MOVIL 290L (75 GAL) - 220V	
G.18951/340	MOVIL 340L (90 GAL) - 220V	
G.18951/380	MOVIL 380L (100 GAL) - 220V	
G.18951/400	MOVIL 400L (105 GAL) - 220V	
G.18951/1000	STATIONARY 1000L (260 GAL) - 220V	
G.18951/1140	STATIONARY 1140L (300 GAL) - 220V	


4.2. Draw and dimensions models working in 220VAC

Models G.18951/A & G.18951/A-1


Models G.18951/B & G.18951/B-1


Models G.18951/C


Model G.18951/D

Model G.18951/E

Model G.18951/X (MADE TO MEASURE)

Bathroom made to measure requested by the client.

4.3. Models working in 120VAC

WHIRLPOOL BATHS 120V	BAÑOS ESTANDAR	
	G.18951/A (120V)	MOVIL 60L (15 GAL) - 120V
G.18951/A-1 (120V)	ADJUSTABLE HEIGHT 60L (15 GAL) - 120V	
G.18951/B (120V)	MOVIL 80L (22 GAL) - 120V	
G.18951/B-1 (120V)	ADJUSTABLE HEIGHT 80L (22 GAL) - 120V	
G.18951/C (120V)	MOVIL 100L (27 GAL) - 120V	
G.18951/D (120V)	MOVIL 170L (45 GAL) - 120V	
G.18951/E (120V)	MOVIL 320L (85 GAL) - 120V	
G.18951/X (120V)	BAÑOS A MEDIDA	
G.18951/40 (120V)	MOVIL 40L (10 GAL) - 120V	
G.18951/230 (120V)	MOVIL 230L (60 GAL) - 120V	
G.18951/290 (120V)	MOVIL 290L (75 GAL) - 120V	
G.18951/340 (120V)	MOVIL 340L (90 GAL) - 120V	
G.18951/380 (120V)	MOVIL 380L (100 GAL) - 120V	
G.18951/400 (120V)	MOVIL 400L (105 GAL) - 120V	
G.18951/1000 (120V)	STATIONARY 1000L (260 GAL) - 120V	
G.18951/1140 (120V)	STATIONARY 1140L (300 GAL) - 120V	

4.4. Draw and dimensions models working in 120VAC

The draws and dimensions of the whirlpools working in 120VAC, are exactly the same than all those working in 120VAC. Please see 4.2.

5. Classification according to Medical Devices Directive 93/42/EEC and normative EN 60601-1:2006

Directive 93/42/EEC

Medical Devices are classified into Class I, IIa, IIb and III. Classification is set out in Annex IX of Directive 93/42/EEC and its amendment 2007/47/EC according to the following general rules.

Rule 9 says:

All active therapeutic devices intended to administer or exchange energy are in Class IIa unless their characteristics are such that they can manage the human body energy exchange with the same potentially hazardous way, taking into account the nature, the density and point of application of energy, in which case they are in class IIb. All active devices intended to control the operation of active therapeutic devices in Class IIb, or intended directly to influence the performance of such devices are in Class IIb.

According to that Annex, the whirlpool baths fall within Rule 9.

En 60601-1:2006 (6)

According to point 6.2 is Class I

According to point 6.3 is IPX0

According to point 6.6 is Working mode= continuous

EN 60601-1:2006 (7.2, 7.9)

According to point 7.2.2. the product is clearly identified in 4.1 and 4.2 of this document. In 14.4 you can find the label that will finally appear in the product. The label is the same than for all models with the only difference of model number and electrical features.

According to 7.2.3 the symbol  when it is necessary in the document

According to point 7.2.6. the product is clearly marked. In 14.4 you can find the label that appears in the product. The label is the same than for all models with the only difference of model number and electrical features. The unit shows a label in the same part of the mains connection with this information (it is shown the label for 120V and 120V):

220 V ~ 50 Hz 3,49A	120 V ~ 60 Hz 7,4A
---------------------------	--------------------------

According to 7.2.10 the symbol  appears in the external surface of product in a label like this

According to 7.9.2.15 it is clearly stated that:

- There is no risk when the product reaches the end of the life cycle. Materials are Inox Steel for the major part of the product. No materials or components could damage to people or goods.
- Ortotecnica always suggest to get in contact with the Wastes Management Entities in your country to dispose correctly the product when it is in the end of the life cycle

According to 7.9.2.16 it is clearly stated that:

- Classifications and markings: please see 6 and 7

6. Operator skills and training

SKILLS

Operators using the whirlpool need:

- ▶ To have knowledge and experience of aquatic physical-therapy procedures.
- ▶ The ability to assist the patient.

TRAINING

Operator trainees need to:

- ▶ Be trained in aquatic-therapy protocols.
- ▶ Be familiar with the types of patients who should or should not receive this type of physical therapy (see Patient Evaluation, below).
- ▶ Read and understand this manual (and the manual of the turbine, if used).
- ▶ Be trained on the use of the whirlpool (and the manual of the turbine, if used).
- ▶ Practice with the whirlpool (and the manual of the turbine, if used) before using with patients.

7. Patient evaluation

Read Working Instructions

Patients who are electrically-susceptible (patients with exposed, non-waterproof electric leads, monitors, etc.), patients carrying infectious disease, or patients with certain other medical, mental or physical conditions should not receive treatment with the whirlpool.

The trained operator must evaluate and verify that each patient is suitable for hydrotherapy treatment before permitting the patient to begin hydrotherapy. If in doubt, consult a medical professional before providing hydrotherapy treatment.

Read Working Instructions

To reduce the risk of exposure to bloodborne diseases such as VIH-1 and hepatitis when using the whirlpool, read and follow disinfecting and cleaning instructions in this manual thoroughly.

8. Inspecting the whirlpool bath

This product has been carefully packaged at the factory to minimize the possibility of damage during shipping.

- ▶ Inspect the packaging for external signs of damage.
- ▶ Inspect the contents for damage.

If there is visible damage to the instrument upon receipt, inform the shipping company and Ortotecsa immediately.

Read Working Instructions

Do not attempt to operate this equipment if there is evidence of shipping damage or you suspect the unit is damaged. Damaged equipment may present additional hazards to you. Contact Ortotecsa technical support for advice before attempting to plug in and operate damaged equipment.

Have your facility's equipment maintenance personnel inspect the whirlpool regularly. Follow the checklist and operate the whirlpool through all its functions as described in this manual.

Inspection Check List

1. Are all components present?
2. Is the whirlpool free of excessive wear?
3. Does the turbine mount securely in place?
4. Is a properly-grounded and voltage-matched hospital grade receptacle available for the turbine?
5. Is the outlet or turbine cord equipped at the desired Height?
6. Does the drain valve open and close properly?
7. Is the thermometer present and is legible?
8. Do the installed accessories operate without interfering with whirlpool bath use or turbine operation?

9. While using the whirlpool

Read Working Instructions

Operators using the whirlpool need

- ▶ Personnel who will work with the whirlpool need to read this manual.
- ▶ Have a plumbing professional install the whirlpool as instructed in installing section. A rough-in drawing is available to assist the installer. Contact Ortopedica Customer Service, bottom of page, for additional installation not covered in this manual. Confirm that the whirlpool operates properly. See Inspecting the Whirlpool section.

General Guidelines for Use

- ▶ Medical advice is beyond the parameters of this manual.
- ▶ The whirlpool is for professional use only.
- ▶ The trained operator must evaluate and verify that the patient is suitable for hydrotherapy treatment before beginning hydrotherapy with whirlpool. See Patient Evaluation section.
- ▶ Follow your hydrotherapy procedures if the patient has a physician's order for treatment (see Patient Evaluation section).
- ▶ Stay with patient at all time.
- ▶ It is the operator's responsibility to ensure safe practices for the patient and themselves.
- ▶ If a turbine is used: Plug the turbine cord only into a receptacle that is voltage matched, properly grounded and polarized. Verify that the receptacle has GFI protection or order a GFCI plug for the turbine cord. Keep patient hair, gown strings and other loose items away from the impeller housing to avoid entanglement and injury. Keep the area around the turbine clear. The turbine requires a minimum of 24" (610mm) of clearance.

Positions of patient and caretaker

- ▶ Depending of the member to be treated, the patient should sit in such way to have a comfortable position. If it is necessary put a towel between body and edge, in case of arm's treatment, for avoiding pain for supporting the arm right on the edge of whirlpool.
- ▶ Caretaker must provide the more comfortable position to the patient. Caretaker is not intended to be in a specific position during treatment.

10. Instructions for use


Read Working Instructions

This product is intended to be used in hydro-therapy procedures. Always has to be used following the indications of this point.

In addition you must to know that

- ▶ The product is intended to be used by people has knowledge and experience of aquatic physical-therapy procedures.
 - ▶ The product is intended to be used by people with the ability to assist the patient.
- The product just has one function and it is provided by activating the element number 9. When it is deactivated the product stops to provide the function.

- | | |
|-----|------------------------------|
| 1. | Brake wheels |
| 2. | Tube of turbine |
| 3. | Guide tube |
| 4. | Knob |
| 5. | Nozzle |
| 6. | Intake |
| 7. | Pressure control valve |
| 8. | Plug |
| 9. | Switch |
| 10. | Drain Valve (Ball valve 1 ") |
| 11. | Thermometer (double scale) |

How to start

1. Determine the location of the whirlpool and then, in the desired location, pressing the brakes wheels (1).
2. Fix the tube of turbine (2) in the guide tube (3), tightening the knob (4) to ensure the desired height.
3. Connect the Nozzle (5) into the intake (6)
4. Turn the pressure control valve (7) fully clockwise, then fully to the left to lubricate the plunger assembly.
5. Connect to the mains pin (plug grounding) (8).
6. Move switch (9) to On position to run the product. For stopping the product, move switch (9) to Off position.

Filling procedure

1. Fill the bath using the fill hose. The water level must be a maximum of 10 cm below the edge of the bath and a minimum of 15 cm above the turbine.
2. Turn the pressure control valve (7) fully clockwise. This provides maximum pressure and agitation. Turn the switch (9). Adjust the pressure of the pressure control valve (7) to achieve the desired level of pressure and agitation.
3. To be drained, turn the drain valve (10) counterclockwise. Drain tube \varnothing 32mm (note, this valve (10) is charged on the side of the bath according to the model of the bathroom. The drain valve can be found in the lower part of the tank depending of the model)
4. If clearing is carried out in a sink, turn the pressure control valve (7), down the turbine to the bottom of the tub, place the nozzle (5) in the sink and turn on the switch (9). To stop the drain switch-off (9). To clean the waste water using the drain valve or sponge / towel.
5. To replace the thermometer (11), loosen the screws located on the back side by rotating in the opposite direction to clockwise to release it. Enter the new (first probe) and secure by tightening the screw again.

11.Cleaning and disinfecting instructions

Cleaning and Disinfecting the Whirlpool

After draining the Whirlpool:

1. Verify that the turbine is OFF.
2. Remove any loose objects or debris that may be present in whirlpool.
3. To Disinfect Ortotecsa recommends that you disinfect the whirlpool after each patient. Spray or wipe all surfaces of the whirlpool tank and thermometer with a hard surface disinfectant, following the disinfectant manufacturer's instructions.
4. Rinse away disinfectant residue with warm, clean water. If needed (or after the last patient of the day) wipe the inside and outside surfaces of the whirlpool with a towel to reduce streaking and spotting.
4. To Clean as needed, or weekly at minimum. Hand-wash all surfaces of the whirlpool and thermometer with a non-abrasive stainless steel cleaner. Scrub only in the direction of the polish lines on the whirlpool. Follow the manufacturer's instructions for use. Rinse with cold water. Wipe inside and outside surfaces with a towel to reduce streaks and spotting.

Disinfecting and Cleaning of the Turbine

1. Drain the tank.
2. Place an empty bucket beneath the turbine. Fill with water and add disinfectant solution per manufacturers recommended use.
3. Turn the control valve to the slow position.
4. Turn on the turbine for 5 seconds and turn it off. Allow the turbine to soak for a 10 minute contact time.
5. Spray or wipe all exposed surfaces with a mixture of suitable disinfectant solution per manufacture's recommended use. Exposed surfaces should remain wet for at least 10 minutes.
6. Drain the turbine bucket into the tank and open drain.
7. Rinse all exposed surfaces with warm/hot water.
8. Fill the bucket with water and run turbine for five seconds to rinse impeller cavity.

Ortotecsa recommends the use of Hydrochlor® Whirlpool Antiseptic with all patient treatments. This simplifies cleaning of equipment by eliminating need for repeated rinsing and scrubbing.

Disinfecting and Cleaning of stainless steel

Clean weekly or more often, as needed (especially high-polish surfaces)

Recommended cleaning materials

Sponge natural or artificial
Nylon or other soft-bristle material brush
Soft cloth (as used on automobile finishes).

Recommended cleaning Solutions

Hand dishwashing liquid/soft water solution
Mild soap/soft water solution
3M Stainless Steel Cleaner/Polish
White vinegar/soft water solution (for brightening, removing oil deposit and hard water deposits)
CLR Brand Cleaner or baking soda /soft water solution (for brightening, removing hard water deposits)
Club soda and sponge

To remove smudges and fingerprints:

Wipe surfaces with a quality Stainless Steel Cleaner / Polish. Apply using a soft non-abrasive cloth. Always follow cleaner/polish product directions provided.

To remove dirt and debris

Wash surfaces with a mild liquid soap. Apply using a soft, non-abrasive cloth. Rinse surfaces thoroughly with clean water. Afterwards, using a soft non-abrasive cloth, wipe surfaces with stainless steel cleaner/polish.

Special Note: After cleaning for serious problems, let dry, and expose to air for at least 24 hours to allow "healing" (restoring of the chromium oxide layer) of the stainless steel surface.

12.Maintenance and troubleshooting

WARNING: DO NOT MODIFY THIS PRODUCT WITHOUT MANUFACTURER PERMISSION

12.1. Maintenance Schedule

The whirlpool requires regular maintenance. Set up and follow a maintenance schedule. The table represents minimum intervals for maintenance.

MINIMUM MAINTENANCE INTERVALS	EACH USE	AS NEEDED	EACH WEEK	EACH MONTH
DISINFECTING	X			
CLEANING		X	X	
INSPECTING		X		X

When using maintenance products, follow the manufacturers' directions and read the manufacturer's MSDS, Material Safety Data Sheets.

12.2. Maintenance procedures

Motor

Motors are equipped with sealed greased ball bearings that do not require lubrication. Vent holes in the motor and motor cover shell must be cleaned periodically to allow free air flow and prevent overheating.

Hard Water

The mineral properties of hard water, if allowed to build up over a period of time, have the potential for creating damage to the moving parts of the turbine assembly. We recommend the installation of a water softener or if that is not possible, the use of a water softening agent.

12.3. Troubleshooting

When there is a significant loss in pressure at the turbine assembly

- ▶ Check water inlet holes underneath the pump, clean away any foreign material.
- ▶ Add water, pump body must be completely submerged.

Motor does not start when switch is turned on

- ▶ Check your power source.

- ▶ Have a qualified technician check switch.
Motor ceases operation
- ▶ Check for motor overloads: Possible causes are:
 - a. Blocked motor ventilation holes
 - b. Excessive current load

Read Working Instructions

The motor is equipped with an overload protector. In an overheated condition, the motor will automatically shut off until it has returned to a normal operating temperature.

To treat that condition properly:

- ▶ Personnel has to have the necessary training
- ▶ It has taken the necessary time to cool the unit
- ▶ The unit has to be disconnected from Mains

12.4. Labelling

Models are described in 4.1 and 4.3.

The labels below are an example for the models G.18951/A(240V) and 18951/A(240V), because all the models has the same label, except the model, that changes according capacity of tank.

Meaning of symbols

Part in contact (whirlpool)

Manufacturer data

Manufacture date (year/month)

SN:

Read instructions before using

Ortecsa:

Serial number

Model:

Trademark

REF:

Model

CE

Reference of the product in the catalogue

CE marking

Electrical features:

120 V/50 Hz/550W

Electrical features:

120 V/60 Hz/550W (models 120VAC)

Earth connection

13. Information related to EMC

EMC EN 60601-1-2-2007

Emissions

The unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

The unit is suitable for use in all establishments, including domestic establishments and those directly connected to the mains supply low voltage which supplies buildings used for domestic purposes. It doesn't built any active electronic component.

Immunity

Floors should be of wood, concrete or ceramic. If floors are covered with synthetic material, the relative humidity should be at least 30%.

The unit is tested to be immune to electric fields of:

- ▶ ± 6 kV by contact.
- ▶ ± 8 kV by air.

The quality of the Mains should be that of a typical commercial or hospital. Normally lamps will support without show any problema of working:

- ▶ Quick bursts: ± 2 kV for Mains.
- ▶ Shockwaves: ± 1 kV differential mode

The quality of the power supply must be that of a typical commercial or a hospital. If the user of the unit requires continued operation during power interruptions, it is recommended that the lamp be powered from an uninterruptible power supply or a battery.

Magnetic fields at Mains frequency should be like a typical location in a commercial environment or hospital: 3 A / m.

Mobile equipment and RF communications should not be used closer to any part of the lamp, including cables. The recommended separation distance to the transmitter frequency is depending of the working freq of the device:

- ▶ $d = 1.17 \sqrt{P}$ 150 kHz to 80 MHz
- ▶ $d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz
- ▶ $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz

Where P is the maximum output power of the transmitter in watts (w) according to the transmitter manufacturer specification and d is the recommended separation distance in meters.

Field strengths from fixed RF transmitter, as determined by an electromagnetic site inspection and should be less than the compliance level in each frequency range of the above table.

Interference may occur in the vicinity of equipment marked with the following symbol: 

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NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation the unit.
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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The whirlpool bath is intended for use in a electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the whirlpool bath can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitter) and the whirlpool bath as recommended below, according to the maximum output power of the communication equipment.
El baño de remolinos está previsto para usarse en un entorno electromagnético en el que se controlan las perturbaciones radiadas de RF. El cliente o usuario del baño puede ayudar a prevenir la interferencia electromagnética manteniendo una distancia mínima entre el equipo portátil y móvil de comunicaciones de FR (transmisores) y el baño según se recomienda debajo, conforme a la máxima potencia de salida del equipo de comunicaciones.

Rated maximum output power of transmitter W Máxima potencia de salida asignada del transmisor	Separation distance according to frequency of transmitter(m)		
	150 kHz to 80 MHz d = 1,17 √P	80 MHz to 800 MHz d =1,17√P	800 MHz to 1 GHz d=2,33√P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,74
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power nor listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 to 80 MHz and 800 MHz, the separation distance is applied in the highest frequency range.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people.

14.Warranty

ORTOTECSA, S.L. will provide a warranty for two years from the date of purchase, the working of this product against manufacturing defects. Our warranty includes repair, replacement, or exchange of the product and / or components including labor in our installations.

To have effect this guarantee, is required the purchase invoice.

THIS WARRANTY SHALL NOT UNDER THE FOLLOWING CONDITIONS:

When the use, care and operation of the product has not been done in accordance with the instructions contained in the operating instructions.

When the product has been used outside of its capacity, abused, beaten, exposed to moisture, ingressed by any liquid or corrosive substance, as well as any other fault attributable to the consumer.

When the product has been disassembled, modified or repaired by unauthorized persons.

The breakage and / or replacement of the bulb are not part of the warranty for this product.

In our continuous improvement process, Ortotecsa reserves the right to modify any characteristic without prior notice.