TANQUE DE HUBBARD

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Última revision 18/06/2019

Contenidos

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			
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation the unit.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Operation and Maintenance manual

Contents

1.	Introduction.....	21
1.1.	About this manual	21
1.2.	Important information	21
2.	Manufacturer.....	22
3.	Purpose of product	22
4.	Model and general dimensions.....	22
5.	Classification according to Medical Devices Directive 93/42/EEC.....	23
6.	Operator skills and training.....	23
7.	Patient evaluation.....	24
8.	Inspecting the Hubbard tank	24
9.	While using the tank	25
10.	Instructions for use.....	25
11.	Cleaning and disinfecting instructions	28
12.	Maintenance and troubleshooting	29
12.1.	Maintenance Schedule	29
12.2.	Maintenance procedures.....	30
12.3.	Troubleshooting	30
12.4.	Labelling	31
13.	Information related to EMC	32
14.	Warranty	37

1. Introduction

1.1. About this manual

This manual is a combined user/operation and maintenance manual.
It is applicable to all models found in point 4.

1.2. Important information

Read Working Instructions

Read these instructions carefully before handling the product.

- ▶ If these instructions are not kept, could result in damage to your new tank and/or body injury.
- ▶ Before connecting the tank, check that your mains voltage corresponds to that is indicated on the label.
- ▶ Do not use the tank if it does not work properly or is damaged, or if the power cable or plug is defective
- ▶ Tank is intended just for professional use only.
- ▶ DO NOT operate tank without properly filling with water.
- ▶ Fill tank with water to the prescribed level before plugging into an electrical receptacle.
- ▶ Operating the appliance without water may result in damage to the motor.
- ▶ Improper use of the tank can cause injury. Use the tank only for the purpose described in this manual.
- ▶ The turbine is top-heavy. Improper handling can cause injury or damage. Handle the turbine with care.
- ▶ Connection improper items to the tank can cause injury to persons and damage to product. Use only Ortotesca approved items in tanks.
- ▶ Never immerse the product in water or use in close proximity to it.
- ▶ Do not use the tank near a bathtub, shower or with wet hands. Should any part of the product was wet, reject the use.
- ▶ For cleaning proceed as it is mentioned in conservation section of this manual.
- ▶ Do not cover the product during use.
- ▶ Don't stack tanks with other products.
- ▶ Use the tank on flat and stable.
- ▶ Keep a safe distance bigger than 80 centimetres, between tank and potentially flammable objects.
- ▶ Turn off the tank and unplug from the mains before cleaning or maintenance.
- ▶ Do not sleep with the product connected.
- ▶ Persons with disabilities, disabled, handicapped or sick children should not use this product without direct adult supervision.
- ▶ It should not be used under the influence of alcohol, sleeping pills or sedatives.
- ▶ If you are carrying a pacemaker or similar cardiovascular device, we recommend you always check with your doctor before using the product.
- ▶ Should you experience pain or unpleasant sensation during use, stop using and consult your doctor.
- ▶ Don't locate the tank for the treatment in such way be difficult unplug it from the outlet if was necessary by fails
- ▶ When the tank has reached its life time should be disposed according to the local government laws. Please get in contact with the local authorities for more information.
- ▶ There are no special environmental conditions for the use. An indoor environment is perfect for the use. Nevertheless the recommendation is: temperature of 25-28°, humidity between 30-50%, Atmospheric pressure: 1000-1025hPa-mb
- ▶ Don't manipulate or open the unit: it is not allowed except by the manufacturer personnel.
- ▶ The Mains cable must NOT be replaced. Only manufacturer is allowed to do it.

2. Manufacturer

Name: ORTOTECSA S.L
Address: C/Peña de Francia,8
28500 Arganda del Rey
Madrid. España
Tlf.: +34 918 719 132
Fax: +34 918 719 323
email: info@ortotecs.com
WEB: www.ortotecs.com


3. Purpose of product

The Hubbard tank is a class of whirlpool, consisting of a tank to immerse the whole body of the patient, whose water is kept in constant agitation by one or two turbines. In these tanks, the thermal effects of hot or cold water are added to derivatives from agitation of stirring water. It is one of the hydrotherapeutic methods most studied and currently used in the treatment of physical dysfunction.

The effects of stirring water on tissues have been studied by several authors. The turmoil created in the tank works as a source of stimulation mechanical skin, acting as a anti-irritant and encouragement of large sensitive afferents, blocking the transmission of pain. Moreover, the increased agitation convective will propagate of hot mechanism.

4. Model and general dimensions

Model	Description
G.18957	TANK 1600 L – 423 GAL
G.18957/A	TANK 1320 L – 349 GAL

5. Classification according to Medical Devices Directive 93/42/EEC

Medical Devices are classified into Class I, IIa, IIb and III. Classification is set out in Annex IX of Directive 93/42/EEC and its amendment 2007/47/EC according to the following general rules.

Rule 9 says:

All active therapeutic devices intended to administer or exchange energy are in Class IIa unless their characteristics are such that they can manage the human body energy exchange with the same potentially hazardous way, taking into account the nature, the density and point of application of energy, in which case they were in class IIb. All active devices intended to control the operation of active therapeutic devices in Class IIb, or intended directly to influence the performance of such devices are in Class IIb.

According to that Annex, the Hubbard tank is Class I by Rule 9.

Classification with EN 60601-1:2006 (6)

According to point 6.2 is Class I

According to point 6.3 is IPX0


According to point 6.6 is Working mode= continuous


Fulfilling with 60601-1:2006 (7.2, 7.9)

According to point 7.2.2. the product is clearly identified in 4.1 and 4.2 of this document. In 12.4 you can find the label that will finally appear in the product. The label is the same than for all models with the only difference of model number and electrical features.

According to 7.2.3 the symbol  when it is necessary in the document

According to point 7.2.6. the product is clearly marked. The symbol  is used in the labels on products. In 12.4 you can find the label that appears in the product. The label is the same than for both models with the only difference of model number and electrical features. The unit shows a label in the same part of the mains connection with this information:

According to 7.2.6 the symbol  is used in the labels on products. See 12.4.

According to 7.2.10 the symbol  appears in the external surface of product in a label like this

According to 7.9.2.15 it is clearly stated that:

- There is no risk when the product reaches the end of the life cycle. Materials are Inox Steel for the major part of the product. No materials or components could damage to people or goods.
- Ortotecsa always suggest to get in contact with the Wastes Management Entities in your country to dispose correctly the product when it is in the end of the life cycle

According to 7.9.2.16 it is clearly stated that:

Classifications and markings: please see 5

6. Operator skills and training

SKILLS

Operators using the tank need:

- ▶ To have knowledge and experience of aquatic physical-therapy procedures.
- ▶ The ability to assist the patient.

TRAINING

Operator trainees need to:

- ▶ Be trained in aquatic-therapy protocols.
- ▶ Be familiar with the types of patients who should or should not receive this type of physical therapy (see Patient Evaluation, below).
- ▶ Read and understand this manual (and the manual of the turbine, if used).
- ▶ Be trained on the use of the tank (and the manual of the turbine, if used).
- ▶ Practice with the tank (and the manual of the turbine, if used) before using with patients.

7. Patient evaluation

Read Working Instructions

Patients who are electrically-susceptible (patients with exposed, non-waterproof electric leads, monitors, etc.), patients carrying infectious disease, or patients with certain other medical, mental or physical conditions should not receive treatment with the tank.

The trained operator must evaluate and verify that each patient is suitable for hydrotherapy treatment before permitting the patient to begin hydrotherapy. If in doubt, consult a medical professional before providing hydrotherapy treatment.

Read Working Instructions

To reduce the risk of exposure to blood borne diseases such as VIH-1 and hepatitis when using the tank, read and follow disinfecting and cleaning instructions in this manual thoroughly.

8. Inspecting the Hubbard tank

This product has been carefully packaged at the factory to minimize the possibility of damage during shipping.

- ▶ Inspect the packaging for external signs of damage.
- ▶ Inspect the contents for damage.

If there is visible damage to the product upon receipt, inform the shipping company and Ortotecsa immediately.

Read Working Instructions

Do not attempt to operate this product if there is evidence of shipping damage or you suspect the product is damaged. Damaged product may present additional hazards to you. Contact Ortotecsa technical support for advice before attempting to plug in and operate.

Have your facilities maintenance personnel inspect the tank regularly. Follow the checklist and operate the tank through all its functions as described in this manual.

Inspection Check List

1. Are all components present?
2. Is the tank free of excessive wear?
3. Does the turbine mount securely in place?
4. Is a properly-grounded and voltage-matched hospital grade receptacle available for the turbine?
5. Is the outlet or turbine cord equipped at the desired height?
6. Does the drain valve open and close properly?
7. Is the thermometer present and is legible?
8. Do the installed accessories operate without interfering with tank use or turbine operation?

9. While using the tank

Operators using the tanks need

- ▶ Personnel who will work with the tank need to read this manual.
- ▶ Have a plumbing professional install the tank as instructed in installing section. A rough-in drawing is available to assist the installer. Contact Ortopedisa Customer Service, bottom of page, for additional information not covered in this manual. Confirm that the tank operates properly.
- ▶ See Inspecting the tank section.

General Guidelines for Use

- ▶ Medical advice is beyond the parameters of this manual.
- ▶ The tank is for professional use only.
- ▶ The trained operator must evaluate and verify that the patient is suitable for hydrotherapy treatment before beginning hydrotherapy with tank. See Patient Evaluation section.
- ▶ Follow your hydrotherapy procedures if the patient has a physician's order for treatment (see Patient Evaluation section).
- ▶ Stay with patient all the time.
- ▶ It is the operator's responsibility to ensure safe practices for the patient and themselves.
- ▶ If a turbine is used: Plug the turbine cord only into a receptacle that is voltage matched, properly grounded and polarized. Keep patient hair, strings and other items away from the impeller in order to avoid injuries.
- ▶ Keep the area around the turbine clear. The turbine requires a minimum of 24" (610mm) of clearance.

Positions of patient and caretaker

- ▶ The whole patient should be inside the tank with the head supported in the head support that it is indicated in the following drawing. The patient should be inside the tank in such way to have a comfortable position.
- ▶ Caretaker must provide the more comfortable position to the patient. Sometimes caretaker must be inside the tank with the patient. The caretaker should adopt those positions that make possible the treatment is going to be provided.

10. Instructions for use



Read Working Instructions

This product is intended to be used in hydro-therapy procedures. Always has to be used following the indications of this point.

In addition you must to know that

- ▶ The product is intended to be used by people has knowledge and experience of aquatic physical-therapy procedures.
- ▶ The product is intended to be used by people with the ability to assist the patient.
The product just has one function and it is provided by activating the element number 9. When it is deactivated the product stops to provide the function.

**How to start**

Locate all the elements before starting

1. Feet anti vibration
2. Tube turbine
3. Guide tube tank
4. Knob
5. Nozzle
6. Intake
7. Pressure control valve
8. Plug
9. Switch
10. Drain Valve
11. Thermometer
12. Thermostatic mixer
13. Fill valve
14. Wrench heat (1 ")
15. Wrench cold (1 ")
16. Sewage input
17. Sewage return

How to install the unit

1. Sure that the tank is level, rectify it it is necessary the leveling legs (1)
2. Fix the turbine tube (2) into the guide tube (3) tank, tightening the knob (4) to ensure the chosen height. To be drained, turn the drain valve (10) counter clockwise. Drain tube Ø 32mm (note, this valve (10) is charged on the side of the bath according to the model of the tank. The drain valve can be found in the lower part of the tank depending of the model
3. Connect the Nozzle (5) in the intake (6).
4. Turn the pressure control valve (7) fully clockwise, then fully to the left to lubricate the plunger assembly.
5. Connecting to the mains (plug grounding) (8).

How to run the unit

1. To run the unit move the switch (9) to On position. To stop working the unit move the switch (9) to Off position.

How to fill the tank

1. Open Wrench heat water (14) and Wrench cold (15).
2. Ensure that the drain valve (10) is closed.
3. Fill the tank using the thermostatic valve (12) and / or the filling valve (13). The water level must be a maximum of 10 cm below the edge of the tank and at least 15 cm above the pump.
4. Turn the pressure control valve (7) fully clockwise. This provides maximum pressure and agitation. Turn the switch (9). Adjust the pressure of the pressure control valve (7) to achieve the desired level of pressure and agitation.

How to drain the tank


1. Verify that the turbine is OFF
2. Remove any loose objects or debris that may be present in tank.
3. Proceed to empty the tank, open the drain valve (10). To clean the waste water using a sponge / towel

11. Cleaning and disinfecting instructions

Cleaning and Disinfecting the Tank

1. To Disinfect Ortoteca recommends that you disinfect the tank after each patient. Spray or wipe all surfaces of the tank and thermometer with a hard surface disinfectant, following the disinfectant manufacturer's instructions.
2. Rinse away disinfectant residue with warm, clean water. If needed (or after the last patient of the day) wipe the inside and outside surfaces of the tank with a towel to reduce streaking and spotting.
3. To Clean as needed, or weekly at minimum. Hand-wash all surfaces of the tank and thermometer with a non-abrasive stainless steel cleaner. Scrub only in the direction of the polish lines on the tank. Follow the manufacturer's instructions for use. Rinse with cold water. Wipe inside and outside surfaces with a towel to reduce streaks and spotting.

Read Working Instructions

Ortoteca recommends the use of Hydrochlor® Tank Antiseptic with all patient treatments. This simplifies cleaning of equipment by eliminating need for repeated rinsing and scrubbing.

Disinfecting and Cleaning of stainless steel

Clean weekly or more often, as needed (especially high-polish surfaces)

Recommended cleaning materials

Sponge natural or artificial
Nylon or other soft-bristle material brush
Soft cloth (as used on automobile finishes).

Recommended cleaning Solutions

Hand dishwashing liquid/soft water solution
Mild soap/soft water solution
3M Stainless Steel Cleaner/Polish
White vinegar/soft water solution (for brightening, removing oil deposit and hard water deposits)
CLR Brand Cleaner or baking soda /soft water solution (for brightening, removing hard water deposits)
Club soda and sponge

To remove smudges and fingerprints:

Wipe surfaces with a quality Stainless Steel Cleaner / Polish. Apply using a soft non-abrasive cloth. Always follow cleaner/polish product directions provided.

To remove dirt and debris

Wash surfaces with a mild liquid soap. Apply using a soft, non-abrasive cloth. Rinse surfaces thoroughly with clean water. Afterwards, using a soft non-abrasive cloth, wipe surfaces with stainless steel cleaner/polish.

To remove rust stains

Wipe surfaces with CRES (available through Whitehall) or equivalent cleaner. Use recommended solutions. Apply using a soft non-abrasive sponge. Rinse surfaces thoroughly immediately after application. Always follow cleaner product directions provided. Afterwards, using a soft, non-abrasive cloth, wipe surfaces with stainless steel cleaner/polish.

For tough problems

CRES Cleaner specifically for rust stains (Available through Ortoteca)

Tarn-X for general stains
7 chrome polish
-Silver polish

To remove stubborn spots or to treat a scratch (Standard Satin Finish only):

Use of synthetic, abrasive, general-purpose pads such as Scotch Brite is recommended. Apply the stainless steel cleaner/polish to the synthetic.

Carefully rub out spot with cleaner/polish. *Be sure to rub in the direction of the grain!* Do not allow steel wool to come in contact with stainless steel. Steel particles can embed into the stainless steel surface and create rust!

Special Note: After cleaning for serious problems, let dry, and expose to air for at least 24 hours to allow "healing" (restoring of the chromium oxide layer) of the stainless steel surface.

For enviro-glaze stainless steel

Protecting the powder coated surface in keeping an Whitehall Enviro-Glazed stainless steel product

For normal wipe-down cleaning, a mild soap and a sponge or cotton cloth will perform well.

Ordinary mild cleaners for home use will sanitize the fixture. (eg. Soft-Scrub, Comet Soft Cleanser Clorox, Scrubbing Bubbles Gel.)

Never use harsh cleansers such as "Kaboom", Lysol Gel, Lysol Cling or "The Works" as these contain powerful acidic or alkaline chemicals which may damage the finish or even the metal of the tank itself. Use of such chemicals will void the warranty.

NEVER use a metal brush while cleaning. The use of such a tool at any time will scratch the finish and perhaps cause corrosion, voiding the warranty.

Things to avoid

Harsh polishing cleanser (Comet, Bon Ami, etc.)

Harsh chemicals (strong acidic or alkaline materials such as hydrochloric acid, sodium hydroxide, etc.)

Any abrasive substance (dirt/mud, ceramic items, etc.)

Chlorinated cleansers, sanitizers or bleach of any kind

Steel wool, steel utensils or any object made of steel or metal, left in contact with stainless surface (hair pins, paper clips, aerosol cans)

Scotch-Brite type pads on high-polished surfaces Seawater or excessively hard water Salt or salty foods allowed to dry on surface

Leaving scratches in metal surface (soil can collect in these, preventing surface from "healing")

Leaving hard water sitting (evaporating) for extended periods, allowing mineral deposits to collect.

12. Maintenance and troubleshooting

WARNING: DO NOT MODIFY THIS PRODUCT WITHOUT MANUFACTURER PERMISSION

12.1. Maintenance Schedule

The tank requires regular maintenance. Set up and follow a maintenance schedule. The table represents minimum intervals for maintenance.

MINIMUM MAINTENANCE INTERVALS	EACH USE	AS NEEDED	EACH WEEK	EACH MONTH
DISINFECTING	X			
CLEANING		X	X	
INSPECTING		X		X

When using maintenance products, follow the manufacturers' directions and read the manufacturer's MSDS, Material Safety Data Sheets. You can purchase a recommended disinfectant from Ortotecs Customer Service (Back of Manual) or your Ortotecs distributor.

12.2. Maintenance procedures

Motor

Motors are equipped with sealed greased ball bearings that do not require lubrication. Vent holes in the motor and motor cover shell must be cleaned periodically to allow free air flow and prevent overheating.

Hard Water

The mineral properties of hard water, if allowed to build up over a period of time, have the potential for creating damage to the moving parts of the turbine assembly. We recommend the installation of a water softener or if that is not possible, the use a water softening agent.

12.3. Troubleshooting

When there is a significant loss in pressure at the turbine assembly

- ▶ Check water inlet holes underneath the pump, clean away any foreign material.
- ▶ Add water, pump body must be completely submerged.

Motor does not start when switch is turned on

- ▶ Check your power source.
- ▶ Have a qualified technician check switch.

Motor ceases operation

- ▶ Check for motor overloads: Possible causes are:
 - a. Blocked motor ventilation holes
 - b. Excessive current load
 - c.

Read Working Instructions

The motor is equipped with an overload protector. In an overheated condition, the motor will automatically shut off until it has returned to a normal operating temperature.

To treat that condition properly:


- ▶ Personnel has to have the necessary training
- ▶ It has taken the necessary time to cool the unit
- ▶ The unit has to be disconnected from Mains

12.4. Labelling

Models are described in 4.1

The labels below are for the models G.18957/A(230V) and 18957/A(115V)

Meaning of symbols

	Part in contact (tank)
	Manufacturer data
	Manufacture date (year/month)
	Read instructions before using
SN:	Serial number
Orteca:	Trademark
Model:	Model
REF:	Reference of the product in the catalogue
CE	CE marking
Electrical features:	230 V/50 Hz/2A
Electrical features:	115 V/50 Hz/2A (models 115V)
	Altern current
	Earth connection
	Class II Product

13. Information related to EMC

EMC EN 60601-1-2-2007

Emissions

The unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

The unit is suitable for use in all establishments, including domestic establishments and those directly connected to the mains supply low voltage which supplies buildings used for domestic purposes. It doesn't built any active electronic component.

Immunity

Floors should be of wood, concrete or ceramic. If floors are covered with synthetic material, the relative humidity should be at least 30%.

The unit is tested to be immune to electric fields of:

- ▶ ± 6 kV by contact.
- ▶ ± 8 kV by air.

The quality of the Mains should be that of a typical commercial or hospital. Normally tanks will support without show any problema of working:

- ▶ Quick bursts: ± 2 kV for Mains.
- ▶ Shockwaves: ± 1 kV differential mode

The quality of the power supply must be that of a typical commercial or a hospital. If the user of the unit requires continued operation during power interruptions, it is recommended that the tank be powered from an uninterruptible power supply or a battery.

Magnetic fields at Mains frequency should be like a typical location in a commercial environment or hospital: 3 A / m.

Mobile equipment and RF communications should not be used closer to any part of the tank, including cables. The recommended separation distance to the transmitter frequency is depending of the working freq of the device:

- ▶ $d = 1.17 \sqrt{P}$ 150 kHz to 80 MHz
- ▶ $d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz
- ▶ $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz

Where P is the maximum output power of the transmitter in watts (w) according to the transmitter manufacturer specification and d is the recommended separation distance in meters.

Field strengths from fixed RF transmitter, as determined by an electromagnetic site inspection and should be less than the compliance level in each frequency range of the above table.

Interference may occur in the vicinity of equipment marked with the following symbol:

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14. Warranty

ORTOTECSA, S.L. will provide a warranty for two years from the date of purchase, the working of this product against manufacturing defects. Our warranty includes repair, replacement, or exchange of the product and / or components including labor in our installations.

To have effect this guarantee, is required the purchase invoice.

THIS WARRANTY SHALL NOT UNDER THE FOLLOWING CONDITIONS:

When the use, care and operation of the product has not been done in accordance with the instructions contained in the operating instructions.

When the product has been used outside of its capacity, abused, beaten, exposed to moisture, ingressed by any liquid or corrosive substance, as well as any other fault attributable to the consumer.

When the product has been disassembled, modified or repaired by unauthorized persons.

The breakage and / or replacement of the bulb are not part of the warranty for this product.

In our continuous improvement process, Ortotecsa reserves the right to modify any characteristic without prior notice.