ELEVADOR / TRANSPORTADOR HIDRÁULICO – G.6500 ELÉCTRICO - G.6505 ELÉCTRICO IP66 – G.6506

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ELECTRICAL – G.6505

ELECTRICAL IP66 – G.6506

Recommended for use in Hospitals, Residential Homes, at home, etc.

Allows for comfortable lifting and lowering of the patient, movement is both effortless and totally secure. Both sides of base open out, using back pedal to allow to place apparatus close to chairs, toilet, or to pass through narrow places.

The apparatus can be easily assembled and put away so that it can be placed in any car boot.

Rotating back casters with brake, (125 mm. diameter) to allow apparatus to be moved more easily when loaded.

Double front casters (75 mm. diameter) to allow for entry into areas raised up to 12 cms. high.

Height of arm adjustable from 180 to 60 cms.

Allowing for patient to be picked up from the floor.

Maximum recommended load: 175 kgs.

The apparatus used in the lift, comply with all the standards required for electrical hoists:

Safety factor 4 in accordance with EN 60601-1. Back-up safety lock. Mechanical clutch release. Maximum stress locking. Emergency release. Noise signal for low batteries.

Supplied with set of batteries, push button control and sling.

IP54 Protection.

ELECTRICAL PATIENT LIFTER IP66 – G.6506

Lift to be in wet areas (swimming pools), Hubbard tank or other types of baths.


DIMENSIONS AND MECHANICAL FEATURES

Dimensions: length: 118cm width: 63cm high: 130cm
Entire lifting arm: 65 to 180 cm
Front wheels: 4 sets of double swivel wheel diam.75mm
Swivel and total lock of diam.125mm: rear wheels.
Opening of adjustable legs by posterior pedal from 62 to 100cm.
Stainless steel frame topped with epoxy.
Recommended maximum load: 175 kgs.

ELEMENTS

- 1 / HAMMOCK
- 2 / HANGER
- 3 / ARM
- 4 / ASAS
- 5 / LIFTING EQUIPMENT
- 6 / MAST
- 7 / PAWS
- 8 / CHASSIS
- 9 / PEDAL
- 10 / BRAKE (REAR WHEEL)


INSTRUCTIONS (Hydraulic Model G.6500)

The arm of the crane will raise or lower powered by hydraulic cylinder.
To vary the spacing of the legs actuate the pedal (9).
For displacements and rotations using the handlebars (4).
To prevent unwanted movement lock the rear wheels (10) by the same brake.

INSTRUCTIONS (G.6505 and G.6506 Electric Models)

It is recommended to set the external charger (14) to the wall by the holder (15) which is attached.
To charge the battery (13) coupling the battery to the charger. The full charge time is about 4 hours (charger has a display indicating battery charge).
To use the lift, attach the battery (13) to the power supply (12) and operate by keypad (16).
The arm of the crane will raise or lower then driven by the actuator cylinder (11).
To vary the spacing of the legs act on the pedal (9).
For displacements and rotations using the handlebars (4).
To prevent unwanted movement lock the rear wheels (10) by the same brake.
IN CASE OF EMERGENCY PRESS THE RED BUTTON of the power supply (12). To unlock that rotate according to the arrows that are on the same button.

LIFTING EQUIPMENT

Lifting equipment installed on electric models (G.6505 and G.6506) consists of the following elements:

11) cylinder actuator 000585:

12) power supply 000605:

Soft Stop.

Audible alarm when its capacity is at 50%.

Electronic overload protection (EOP)

Emergency Button

13) battery 000588

14) external charger 000607

Input voltage: 100 - 240V (50.60Hz)

Current load: up to 650mA

Charging time: 4 hours.


Yellow LED indication of charging

Green LED indicating network connection.

15) support to fix the charger to the wall 000610

16) Remote control 000608

ASSEMBLY

HARNESSES, METHODS OF PLACEMENT:

REMEMBER: BANDS LONG ON LEGS, SMALL SHOULDER

1. LATERAL DECUBITUS

Turn the patient in lateral decubitus (in patients with hip fracture, place a pillow between your legs to prevent gathery), place the harness as in Figure 1 and then turn the patient to the opposite supine position and finish to stretch the harness .

Then the leg bands passed under both legs and intersect, the right leg to hang on the left hook and left leg on the right hook (as in Figure 2)


Figure 1

IT IS VERY IMPORTANT that in patients with hip fracture, the legs stay separated: that's why harness straps intersect. Figure 2.



Figure 2

2. FROM SEATED POSITION

Place the patient harness as in Figure 3. From top to bottom. Crisscrossing straps of the legs as in the previous step

Figure 3

After installing the harness, Raise a bit the arm of the lift close to the bed or chair where is the patient. Place the shorter bands (shoulders) in the back of the hanger, and longer (legs) in the hanger hooks

Gently raise to a height of separation between the bed and the patient, **OPEN THE LEGS OF THE LIFT TO INCREASE ITS STABILITY** and ensure greater security in the movement; place the patient to the chair or vice versa; located him in the chair and lower down to the sitting position; one of the assistants will make sure that the patient is firmly seated pulling the hangers of the harness itself. In the event that the patient may need additional fastening and with the harness on, the longer bands of the legs will be used to secure in the arms of the chair so the patient will not slip out of the chair so easily.

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