TELAR DE MESA COMPLETO -T.O.100

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
COMPLETE TABLE LOOM – T.O.100

The rotation of any of the controls of the warp, makes all the others rotate in a synchronised manner, this synchronisation allows therefore to carry out combined active assisted exercises with both limbs.

For example: flexion and extension of the left wrist, automatically mobilises the right one by flexing and extending, or any elbow, forearm or wrist movement or pronation/supination.

This loom has the appropriate mechanism for placing resistance and is synchronised as mentioned in the previous movements, so that they can be undertaken gradually and progressively, or kept constant.

Picking up the inlaid yarn is carried out by the horizontal displacement of the approach and separation of the fulling device (the frame of the comb), a movement that facilitates active or active assisted flexion and extension of the elbow and shoulder.

The freeing and locking of the previous and following reels, that carry the weave or the loom takes place using the front and side controls respectively.

OPTIONAL ACCESSORIES:

WARP PREPARATION FRAME – T.O.103-1


LOOM THREAD KNIFE – T.O.103-3

LOOM SHUTTLE NEEDLE (30cms.) – T.O.103/A

LOOM SHUTTLE NEEDLE (40cms.) – T.O.103/B

LOOM SHUTTLE NEEDLE (50cms.) – T.O.103/C

LOOM SHUTTLE NEEDLE (60cms.) – T.O.103/D

USER MANUAL

USER MANUAL

1. Place the loom on a solid table.
2. Loose the breaks F1 and F2 to release the reel-plate thread holder C1 and C2 anterior and posterior.
3. Put the threads in the rear plate C2 passing them through the comb (P) and the needle (L) and fixing them on the plate of the previous reel C1.
4. Tight the threads by turning the knobs on the reel s C1 and C2, fixing its position by tightening the brakes F1 and F2.
5. The prono-supination can be made:
 - a. Rotating the wooden handles E1 and E2 in right/left direction
 - b. Rotating the flywheels V1 and V2
 - c. Rotating the rollers R1 and R2
6. Introduce the shuttle needle right-left, perform the prono-supination and reinsert the shuttle needle left-right; repeat the prono-supination to weave.
7. Move the comb P1 and P2 in the direction of the arrow to compact the weave.
The pressure of the flywheels and the rollers can be graduated by turning the knobs G1 and G2.