# CALENTADOR DE PARAFINA (22L.)

## G.18907

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# PARAFINE HEATER (22L.) G.18907

## Revision History:

Rev	Reason	Date
1	Creation of the document	17- jun - 2013
5	Drainage recommendation	04- May - 2018

<b>INDEX</b>	<b>PAGE</b>
MODELS	14
EQUIPMENT CHARACTERISTICS	14
ELECTRICAL SPECIFICATIONS	14
WARNING	16
INDICATIONS OF PARAFFIN	17
CONTRAINdicATIONS AND PRECAUTIONS	17
ELEMENTS	19
INSTRUCTION FOR USE	19
APLICACION	20
PRODUCT STORING	20
CLEANING PRODUCT	21
TROUBLESHOOTING GUIDE	21
LABEL ON PRODUCT	22
WARRANTY	23

## Models

<u>REF</u>	<u>NAME</u>
<b>G.18907</b>	Paraffine bath for Mains Input Voltage of 230Vac

## Equipment characteristics

Paraffin bath with heating system based on "Bain Marie". For hospital use.

Mobile (lock wheels).

Built both the inner tank and the outer bowl in stainless steel AISI 316.

Edge insulation (wood).

Safety mechanism against boiling based on Thermostat.

Thermostat from 0 to 60 °. Thermostatic control paraffin maintains the optimum temperature for maximum therapeutic benefit.

Front panel provided with pilot and thermostat control.

Optional accessories:

- G.18907 - 1: paraffin bath trolley
- G.18907 - 2: filter for paraffin
- G.18910: 1kg paraffin tablet

## Electrical specifications

### Rating

Models		G.18907
Voltage	110 Vac ≈	230 Vac ≈
Frequency	60 Hz	50 Hz
Power	1800 W	1800 W
Amperage	18A	9,1A
Useful life	10 Years	10 Years

### Main cordon:

- 4 m longitude (13 ft)
- Cordon and plug with ground connection.

## Environment Conditions

### Functional environments conditions recommended

- Temperature: 10 °C - 35 °C (Temperature recommended: 20 °C - 25 °C)
- Humidity: 35 - 75%
- Pressure: 1000 - 700 mBar

### Storage and transport Conditions

- Temperature: -5 °C - 70 °C
- Humidity: 0% - 90%

- Pressure: 1000 - 100 mBar

**Environment requirements**

- according to directive 2011/65/UE (RoHS2).
- according to directive 2002/96/EC. (WEE).

**WARNING!**

Read the instructions carefully before handling the product.

- Before connecting the product, check that your mains voltage corresponds to that is indicated on the label.
- Do not use the equipment if it does not work properly or is damaged, or if the power cable or plug is defective.
- Never immerse the product in water or use in close proximity to it.
- Do not use the equipment near a bathtub, shower or with wet hands. Should any part of the product was wet, reject the use.
- For cleaning proceed as it is mentioned in conservation section of this manual.
- Do not cover the equipment during use.
- Dont stack bath with other equipment.
- Use the bath on flat and stable, preventing the unit from falling.
- Keep a safe distance bigger than 80 centimeters, between baths and potentially flammable objects.
- Turn off the bath and unplug from the mains before cleaning or maintenance.
- No treatment must applied after taking analgesics or painkillers, ther heat sensitivity may be diminished by the effect of drugs and this can cause skin lesions.
- The ingestión of alcohol or drugs can also alter the sensitivity to heat.
- The device is not intended for persons insensitive to heat.
- Excessive exposure to heat can cause burns or skin lesions.
- Do not use the paraffin in areas of the body with burns, infections, inflammations, rashes, wounds or hypersensitive areas.
- Do not sleep with the product connected.
- Persons with disabilities, disabled, handicapped, insensitive to heat, or sick children should not use this product without direct adult supervision.
- It should not be used under the influence of alcohol, sleeping pills or sedatives.
- If you are carrying a pacemaker or similar cardiovascular device, we recommend you always check with your doctor before using the product.
- Should you experience pain or unpleasant sensation during use, stop using and consult your doctor.

**INDICATIONS OF PARAFFIN:**

The wax used in therapy must be white, odorless, tasteless and solid and is supplied in the form of plates. Use Medical paraffine.

Therapeutic effects are found with the application of heat:

- To Increase connective tissue extensibility.
- To decrease joint stiffness.
- Analgesic effects
- In subacute and chronic inflammatory processes. Example, rheumatoid arthritis, bursitis, tendinitis and periarthritis escapulahumeral, because it provides a long antiinflammatory and analgesic action.
- Muscle relaxant effect.
- In painful processes as sacrolumbalgias, back pain, cervical pain, sciatic pain, herniated discs, etc..
- Antispasmodic effects.
- Anti-inflammatory effect.
- The efficacy of thermotherapy applications is clinically documented to fight against muscle spasms of defense, subacute and chronic traumatic events and other processes like sympathetic dystrophy - reflects, Raynaud's disease, pelvic inflammatory disease, etc..

**CONTRAINDICATIONS AND PRECAUTIONS****Disorders in the sentitivity.**

Before applying, you must evaluate the thermal and pain sensitivity of the area as well as the patient's circulatory status, because the heat is generally contraindicated or should be used with particular caution on areas anesthetized or obtunded patient.

In most applications, dosimetry is not exact, and must be taken of the patient perceives discomfort or pain when warming threshold was exceeded.

**Severe circulatory disorders.**

Heating of tissues with inadequate blood supply is contraindicated, as the temperature elevation increases the metabolic demands, without associated vascular adaptations, with the risk of producing an ischemic necrosis.

Circulatory disorders, especially arterial, can provide to a decrease in convective mechanism supplied heat dissipation, with the risk of burning the tissues of the area. Altered vessels may have an increased contractility and their response to intense thermal stimulus is essentially spasmodic.

**Bleeding disorders.**

Dont treat areas where there is a tendency to bleed, because the thermal stimulus produces an increase in blood flow. This is the case in diseases such as hemophilia, in post-traumatic states or in patients undergoing long-term steroid therapy, in which there is a tendency to capillary fragility.

**Malignant neoplas.**

Except for special treatment (hyperthermia anticancer), thermotherapy is contraindicated in areas where malignancy is localized, as subtherapeutic temperatures, as used in physiotherapy, can increase the rate of tumor growth. There are many studies showing that thermotherapy applications may increase the development of metastases.

**Pregnant.****Febrile states.**


Avoid applications over recent injuries, whether or not infected by the production risk of burns. Be very careful in the application of recent and skin grafts, and, similarly, on infected skin lesions (risk of exacerbation of these) in peripheral vascular disease or any condition that affects the distribution of internal heat.

It is important to supervise  
the level of the water and that  
it is the proper one.

**ELEMENTS:**

1. Tank of Paraffin
2. Thermostat
3. Pilot
4. Switch to ON/OFF
5. Water level indicator
6. Wheel brake
7. Mains cable connection
8. Hole to fill the water tank
9. Cap water tank
10. Lid of paraffin tank
11. paraffin tank drain
12. drain water tank
13. Reset button

**INSTRUCTION FOR USE:**

Before connecting the unit to the network is necessary to fill the tank with water.

Fill the water tank (8). Check that the tank is full to capacity by the water level indicator (5), it must be filled completely.

With the water tank full of water, fill the tank (1) with paraffin, insert the plug (7) to the Mains and turn the switch (4) to the position 1 (On).

The pilot light (3) will ON.

Select the temperature for the paraffin desired for the thermostat (2).

When the pilot turns OFF, it had reached the temperature selected.


We recommend a 50-55°C temperature.

The warning pilot, will blink from time to time, indicating that the resistance that heats the water in the tank is working to maintain the temperature selected on the thermostat.

Both models of bath are equipped with a system that cuts power before the water reaches 100 ° C. in case of failure of the main thermostat.

If Mains accidentally was cut, the unit can be reset after a period of about 15 minutes for it to access the bottom of the unit, and press the reset button (13).

If a sign of failure arises again, in the main thermostat, so you should call Customer Service.

**Attention!**

Do not stack to another equipment or device.


Using distilled water for good maintenance of the paraffin bath is recommended.

To proceed to the emptying tank, access outfalls paraffin (11) and water (12).

**Caution:** Periodically check the water level. The implementation team with insufficient water level can damage the heating elements.

SELECTOR FOR TEMPERATURE  
OF PARAFFINE. Thermostat.

**CONTROL PANEL**

MAIN SWITCH

 PILOT FOR INDICATING TEMPERATURE  
HAS BEEN REACHED

**Important Notice!:** The temperature of the wax will be selected by the thermostat (2), does not correspond with the temperature dial indicator water level (5), since it does not act as a thermometer (this will mark a difference to least about 20 °).

**APPLICATION**

- Before treatment, the body part should be cleaned with soap and water, then with alcohol to remove any soap residue and prevent bacterial growth on the bottom of the bath.
- Carefully pour the wax with a casserole on the hand or foot. Then take off, in order to form a thin coating of paraffin slightly hardened and adherent to the skin .
- The operation is repeated 8 to 12 times until it forms a thick layer of paraffin wax. Then the area is wrapped in a plastic bag and covered with several towels in order to facilitate heat retention.
- The body area should be removed of any metallic objects and must ensure that the area does not move once introduced in the bath to avoid the appearance of "hot spots".
- By moving the part introduced in the molten paraffin, there is a risk of interrupting the semisolid paraffine barrier, by which the patient could feel a burning sensation. The hand will be immersed with fingers extended and separated as possible.
- The patient should be in a comfortable position, with elevated until the end of the treatment, to avoid the potential appearance of edema.
- The application is maintained for 15 to 20 min.
- After this time, remove the towel and bag, and remove the layer of solid paraffin and put into the bath, but previously sterilized at 80 ° C, or with boiling water
- After application, verify the condition of the skin.
- After application, the area should be cleaned with soap and water. The cleaning can be completed with a gentle massage with a moisturizing lotion or mineral oil to moisturize and soften the skin.
- After applying wax, the skin is smoother, softer, moist and flexible, making it easier to massage and mobilization.

**Important notice:** This product is not intended to submerge any part of the body inside it. In such a case, Ortotecsa is not responsible for its improper use.

**PRODUCT STORING**

- To extend the life of the product, you must unplug it and keep it disconnected from the outlet when not be used for some periods of time.
- Store in a dry place.
- Allow to cool before storing.

- Don't stack objects on top.
- Do not place heavy objects.

**CLEANING PRODUCT**

The bath needs continued supervision to avoid contamination. The bath should be cleaned regularly, as recommended by ORTOTEC莎. The times you reuse the paraffin determine the frequency with which the bath is to be cleaned, but it is recommended in intervals not exceeding six months.

**MAINTENANCE GUIDE**

1. Inspect the environmental conditions in which the bath is located. Avoid storage in wet places. Do not stack other devices or equipment on top of the paraffin bath.
2. Perform integral external and internal cleaning service.
3. Externally inspect the equipment and its main accessories.
4. Check cord and other electrical components.
5. Lubricate moving parts with any commercial lubricant spray.
6. Check the accuracy of temperature control.
7. Check the correct operation of the equipment in conjunction with the operator.
8. The glove or sock which is formed during the treatment should not be added directly to the bath, but collect and sterilized at 80 ° C, or with boiling water.
9. It is also necessary to renew entirely the contents of the bath every six months and clean it properly, since the surface with traces of skin is an excellent medium for bacterial culture.
10. In physiotherapy services it is recommended to have at least one unit for hands and one for feet.
11. The unavoidable dirt residues originate by the paraffin in the surrounding zones and in table, on the floor etc., should be put in a separate area.

**PARAFFIN CLEANING SYSTEM**

- Paraffin cleaning is performed according to their use and degree of dirtiness.
- The temperature is increased to proceed to the emptying of the bath (Recommended 90°).
- Empty the bath content in a tank for this type of residues according to the laws in force in the country. Don't pour down into the drain or down the toilet.
- The procedure to clean the bathroom with soap and water. It should be necessary to pay attention in all those areas with corners and more difficult to reach. If it is necessary use a brush.

**TROUBLESHOOTING GUIDE:**

Possible malfunctions.

INCIDENCE	RECOMMENDATION
Bath doesn't switch on	<p>Check that there is right voltage at the socket used. Press the switch to switch-ON .</p> <p>Check that the thermal cut-out has not tripped (this can be caused by paraffin overheating or caused for not maintaining the water level of the tank at the level indicated in the instructions).</p> <p>Check that the equipment is used with the proper water level (otherwise, it can trip the thermal cut-out and burn resistance). Check the conditions of the power cord.</p>

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**LABEL ON PRODUCT**

**MEANING OF SYMBOLS**

SYMBOL	DESCRIPTION
	MANUFACTURER
	MANUFACTURE DATE (YEAR/MONTH)
<b>REF:</b>	REFERENCE MODEL OF DEVICE
<b>LOTE:</b>	PART NUMBER OF DEVICE AND SERIAL
	CONSULT THE ACCOMPANIMENT DOCUMENTATION
	The device must not be disposed as an urban waste, because it is an Electrical and electronic waste

<b>WARRANTY</b>
-----------------

ORTOTECSA, S.L. will provide a warranty for two years from the date of purchase against manufacturing defects. Our warranty includes repair, replacement, or exchange of the product and / or components including labor in our installations.

To have effect this guarantee, is required the purchase invoice.

To have effect this guarantee, is required the purchase invoice.

**THIS WARRANTY SHALL NOT UNDER THE FOLLOWING CONDITIONS:**

- When the use, care and operation of the product has not been done in accordance with the instructions contained in the operating instructions.
- When the product has been used outside of its capacity, abused, beaten, exposed to moisture, ingressed by any liquid or corrosive substance, as well as any other fault attributable to the consumer.
- When the product has been disassembled, modified or repaired by unauthorized persons.
- The breakage and / or replacement of the bulb are not part of the warranty for this product.

In our continuous improvement process, Ortotecsa reserves the right to modify any characteristic without prior notice.

# ORTOTECSA,S.L.

**Fabricante/Manufacturer: Ortotecsa,s.l.  
C/Peña de Francia,8 – 28500 Arganda del Rey (Madrid), España.  
Tlf.: +34918719132 / Fax: +34918719323  
email: [info@ortotecsa.com](mailto:info@ortotecsa.com), [www.ortotecsa.com](http://www.ortotecsa.com)**