

ESCALERA DE DEDOS / FINGER LADDER

F.42/42

De construcción metálica, rematada en epoxy.
Made for metal with epoxy finish.

F.42/43

Construida en madera estufada y barnizada.
Made from wood treated in the stoving process and varnished.



Manual de uso / User manual

RECOMENDACIONES DE USO / USE RECOMMENDED:

Este dispositivo consiste en una serie de muescas o salientes en los que el paciente debe ir colocando sus dedos. Se utiliza primordialmente para aumentar la amplitud de movimientos en el hombro, sobre todo pidiendo la flexión de miembro superior con el codo en extensión y después llevando a cabo la abducción. Han de realizarse pocas repeticiones pues es un ejercicio muy útil para aumentar la movilidad del hombro, pero excesivamente fatigoso desde el punto de vista muscular.

Siga cuidadosamente todas las instrucciones dadas por su médico, fisio o terapeuta para el uso del producto.

This device consists of a series of grooves or ridges where the patient should be placing his fingers. It is primarily used to increase the range of motion in the shoulder, especially asking flexing the upper limb with the elbow extended and then carrying out the abduction. Few repetitions should carry out as it is a useful exercise to increase shoulder mobility, but too tiring from the muscular point of view.

Follow carefully all instructions given by your doctor, physic or therapist for the use of the product.

RECOMENDACIONES DE INSTALACIÓN / INSTALLATION:

