

PRONOSUPINADOR / PRONOSUPINATOR

SENCILLO - F.42/47-S


Asidero giratorio, provisto de mecanismo opositor de resistencia instalado sobre placa de fijación a pared.

De construcción metálica rematado en epoxy.

STANDARD – F.42/47-S

Rotating handle, equiped with a mechanism to place resistance, fixed on a plaque to be put on the wall.

Metal frame with epoxy finish.



RECOMENDACIONES DE USO:

Consiste en que el paciente agarre el asidero con la mano y la mueva continuamente realizando así la pronosupinación de las muñecas.

Con esta actividad se pretende mejorar la potencia y la resistencia muscular. Si objetivo es rehabilitar, mantener o aumentar la flexibilidad de la muñeca y del antebrazo y su movilidad.

RECOMMENDED USE:

Consists that the patient holds the grip with his hand and move continuously performing pronosupination of the wrist.

This activity is to improve power and muscular endurance. Its objective is to restore, maintain or increase the flexibility of the wrist and forearm and its mobility.

RECOMENDACIONES INSTALACIÓN / INSTALLATION RECOMMENDED:

