

MUELLE PARA CIRCUNDUCCIÓN DE MUÑECA - F.42/48

Asidero de madera barnizada.
Varilla guía cromada.
Bastidor metálico rematado en poliéster.

SPRING FOR THE CIRCUMDUCTION OF THE WRIST – F.42/48

Varnished wooden handle.
Chrome spring rod.
Metal frame with epoxy finish.



INSTRUCCIONES DE USO / USER MANUAL

F.42/48



RECOMENDACIONES DE USO

El **Muelle de Reader** consiste en un circuito con asidero de manos para ejercicios de extremidades superiores, circunducción de muñeca y movimiento del hombro.

Colóquese frente al muelle y cogiendo la empuñadura ir subiendo a lo largo del mismo.

Siga cuidadosamente todas las instrucciones dadas por su médico, fisio o terapeuta para el uso del producto.



RECOMMENDATIONS OF USE

The Reader spring consists of a circuit with hand grip for exercises of the upper extremities, wrist circumduction and shoulder movement.

Face the spring and catching the handle go to rise along the same.

Follow carefully all instructions given by your doctor, physic or therapist for the use of the product.



RECOMENDACIONES DE INSTALACIÓN / INSTALLATION RECOMMENDED

