

ESPEJOS:

MÓVIL - F.90-1

MÓVIL CON ASIDEROS - F.90

MÓVIL CUADRICULADO - F.90-3

MÓVIL CUADRICULADO CON ASIDEROS - F.90-2

ESPEJOS INFANTILES CON ASIDEROS:

MÓVIL - F.91N

MÓVIL CUADRICULADO - F.91N-2

Espejo basculante sobre bastidor metálico. Estructura metálica rematada en epoxy
 Dispone de asideros frontales que sirven de apoyo para permitir al paciente la marcha con el espejo a modo de andador. (F.90 + F.90-2+F.91N + F.91N-2).

Las cuadrículas están grabadas por sistema "chorro de arena". (F.90-2 + F.90-3+F.91N-2)

INSTRUCCIONES:


Montaje:

Colocar el espejo (1) en su soporte (2) intercalando las arandelas de goma (3).

Colocar los pomos (5) intercalando la arandela de acero (4).

Empleo:

Para regular la inclinación, aflojar ligeramente los pomos (5), regular y volver a apretar los pomos.



MIRRORS:

MOBILE MIRROR – F.90-1

MOBILE MIRROR WITH HANDLES – F.90

MOBILE MIRROR DIVIDED INTO SQUARES – F.90-3

MIRROR DIVIDED INTO SQUARES WITH HANDLES– F.90-2

MIRRORS WITH HANDLES INFANTILE:

MOBILE– F.91N

MOBILE WITH SQUARES– F.91N-2

Tilting mirror in a metal frame. Metal frame with epoxy finish.

The square glass models have been cut using the "sand blast" method. (F.90-2 + F.90-3+F.91N-2)

The frame have front handles to help the patient to use the mirror as they walk. (F.90 + F.90-2+F.91N + F.91N-2).

INSTRUCCIONES:

Mounting:

Place the mirror (1) in its support (2) by inserting the rubber packing rings (3).

Place the knobs (5) inserting the steel packing ring (4).

User:

To adjust the angle, slightly loosen the pommels (5), regular and retighten the pommels.

